

# BD™ Blunt Fill and BD™ Blunt Filter Needles

## Lower cost, reduced risk.

BD Blunt Fill and BD Blunt Filter Needles address a wide range of your medication preparation needs. A complement to the full line of BD products, the BD Blunt Fill Needle family provides a substitute for a conventional needle.

The BD Blunt Needle family provides:

- SKU reduction in conventional needles used for preparation
- A lower cost alternative than sharps safety-engineered devices, while removing the need for conventional needles
- Regular and filter needle versions



**18 G 1"**



**18 G 1 1/2"**



**19 G 1"**



**19 G 1 1/2"**



**20 G 1"**



**20 G 1 1/2"**



**21 G 1"**



**21 G 1 1/2"**

BD™ Blunt Fill Needle Family		
Catalog #	Description	Box/Case
<b>BD™ Blunt Fill Needle</b>		
305180	18 G 1 1/2"	100/Box 1000/Case
<b>BD™ Blunt Fill Needle, with BD Luer-Lok™ Tip</b>		
305060	18 G 1 1/2" 3 mL Syringe	100/Box 800/Case
305062	18 G 1 1/2" 5 mL Syringe	100/Box 400/Case
305064	18 G 1 1/2" 10 mL Syringe	100/Box 400/Case
<b>BD™ Blunt Filter Needle</b>		
305211	18 G 1 1/2"	100/Box 1000/Case

**Needle Without Filter**  
(Enlarged to show detail)



**18 G 1 1/2"**

**Easily recognizable color**



**5 Micron Filter**

**Filter Needle**  
(Enlarged to show detail)