

BD™ Oral Syringe System

Recommended for oral medications for patient safety.

BD offers a broad range of syringes and accessories for oral medication administration. Both clear and amber syringes feature mL and teaspoon scale markings for accuracy. The amber syringes meet USP requirements for photo-sensitive medications.²



BD™ Oral Syringe System

Catalog #	Description
-----------	-------------

Clear with Tip Cap

305217	1 mL Syringe
305220	3 mL Syringe
305218	5 mL Syringe
305219	10 mL Syringe

Amber with Tip Cap

305207	1 mL Syringe
305210	3 mL Syringe
305208	5 mL Syringe
305209	10 mL Syringe

Accessories

305222	Bottle Adapter
305223	Oral Syringe Filling Connector

BD Tamper-Tuf™ Medication Containers

BD Tamper-Tuf containers are offered in two sizes (15 mL or 30 mL) providing extra convenience with a large filling port, snap top and volume markings. The containers are available unassembled for non-liquid and liquid oral medication or assembled for liquid medication.

BD Tamper-Tuf™ Medication Containers

Catalog #	Description
-----------	-------------

Assembled

305203	15 mL Container
305204	30 mL Container

Unassembled

305205	15 mL Container
305206	30 mL Container
305228	Container Closer



² Data on file at BD.