

# BD Anesthesia Products

Because Your Patients Trust You to Choose the Best



**FEATURING**  
**Spring Wound Catheters**



Helping all people  
live healthy lives

## BD Epidural and Combined Spinal/Epidural Procedural Trays

BD offers some of the most advanced and clinically preferred technologies available today in Epidural and Combined Spinal Epidural Anesthesia Trays.\*

As a leader in the manufacture of anesthesia needles, syringes, catheters, and other components, BD gives you all the recognized clinical benefits of Spring Wound Catheters... plus a significant savings over other trays.

Choose BD for its total commitment to regional anesthesia through its products and services... because your patients trust you to choose the best.

\*Data on file.

**BD Medical**  
1 Becton Drive  
Franklin Lakes, NJ 07417  
[www.bd.com/anesthesia](http://www.bd.com/anesthesia)

# Spring Wound Epidural Catheters

When you need to balance precision with patient comfort.

The Spring Wound Epidural Catheter enhances patient comfort and safety.\*

- Choice of soft, rounded, closed tip with multiple ports, or open tip available in a variety of tray configurations
- Large, 1" proximal and 2" open series of coils at the distal tip to allow visualization of fluids
- Blended nylon over stainless steel wound spring construction with circumferential centimeter markings



## New Enhanced BD Continuous Epidural and Combined Spinal/Epidural Trays Featuring the Spring Wound Catheter

Component	Continuous Epidural Tray Catalog #									CSE Cat. #
	406016	406039	406113	406116	406117	406118	406038	406114	406115	
25G x 4 11/16 in. BD™ Whitacre Spinal Needle, TW										X
17G x 3.5 in. BD Perisafe™ Weiss Epidural Needle (Fixed Wings)		X	X	X	X	X	X	X	X	X
17G x 3.5 in. BD Perisafe™ Tuohy Epidural Needle	X									
19G x 36 in. Spring-Wound Epidural Catheter OE							X	X	X	
19G x 36 in. Spring-Wound Epidural Catheter CE	X	X	X	X	X	X				X
7 mL BD Epilor™ Syringe LOR Plastic LS	X					X			X	
5 mL BD Epilor™ Syringe LOR Glass LS, Glass tip					X	X	X	X		
5 mL BD Epilor™ Syringe LOR Glass LS, Metal tip		X								X
10 mL BD Epilor™ Syringe LOR Glass BD Luer-Lok™, Metal Tip										
Lidocaine HCL (1.5%) with Epi (1:200,000), 5 mL	X	X	X	X	X	X	X	X	X	X
Lidocaine HCL (1%), 5 mL	X	X	X	X	X	X	X	X	X	X
Sodium Chloride (0.9%), 10 mL	X	X	X	X	X	X	X	X	X	X
18G x 1.5 in. Filter Needle (5 micron), Blunt Tip	X	X	X	X	X	X	X	X	X	X
18G x 1.5 in. Needle, Short Bevel	X	X	X				X	X		X
22G x 1.5 in. Needle	X	X	X	X	X	X	X	X	X	X
25G x 1 in. Needle										X
25G x 1.5 in. Needle			X		X	X	X	X	X	
30G x 1 in. Needle	X									
3 mL Syringe BD Luer-Lok™	X	X	X	X	X	X	X	X	X	X
5 mL Syringe BD Luer-Lok™	X									X
3 mL Syringe BD Luer-Lok™ W/25G x 1" Needle					X	X			X	
10 mL Syringe Luer Slip W/O Tip Protector Latex Free			X		X					
10 mL Syringe BD Luer-Lok™	X							X		
20 mL Syringe BD Luer-Lok™		X						X		X
20 mL Syringe Plastic – LS			X	X		X	X		X	
Sharps Stick Pad	X	X		X			X	X		
Drape, Fenestrated, Blue Paper (4 in. Round, 30 x 26 in.)							X			
Drape, Fenestrated, Clear Plastic (4 in. Round, 32 x 24 in.)	X	X	X	X	X	X		X	X	X

All Spring Wound Catheter Trays contain: BD Perisafe™ Epidural Catheter Connector, 0.22 Micron Flat Filter, Catheter Support Pad, 3 x 3 in. Gauze Pads, Povidone-Iodine Topical Solution, 3/4 oz., Paper Towel, Underpad (19 x 16 in., 3-Ply)

\*Data on file, Third Part Published data on Spring Wound Catheters.  
BD and BD Logo are trademarks of Becton, Dickinson and Company. ©2008 BD.