



BD™ Nerve Block Trays

Because your patients trust you to choose the best



Helping all people
live healthy lives

Premium Trays

- Standard and custom trays designed to meet your procedural needs in the hospital and pain clinic

Precision Components

- BD™ Quincke or BD™ Whitacre Needles available in a variety of lengths and gauges
- Chiba Needles with short bevel point geometry offered in 6" and 8" lengths
- Short bevel and hypodermic bevel needles also available in a full range of sizes

Accessories

- Wide range of accessory components available, including extension tubing sets, control syringes, hypodermic needles, and latex-free syringes

BD™ Nerve Block Trays

Because your patients trust you to choose the best

All Trays Contain:

- Prep section including towel, prep cup, prep swabs, and gauze pads
- Povidone iodine topical solution, USP (1% available iodine)
- Fenestrated drape
- Sharps stick pad (not included in BD single-deck configuration)
- Adhesive bandage (not included in BD single-deck configuration)



Hypodermic Needle



Short Bevel Needle



Block Bevel



Quincke Bevel



Chiba Needle

Material Number	Code Number	20 G x 3.5 in. BD Perisafe™ Tuohy Needle	30 G x 0.5 in. BD™ Hypodermic Needle	27 G x 1.25 in. BD™ Hypodermic Needle	25 G x 1.5 in. BD™ Hypodermic Needle	25 G x 1.0 in. BD™ Hypodermic Needle	18 G x 1.5 in. BD™ Hypodermic Needle	18 G x 1.5 in. BD™ Filter Needle (3 micron), Blunt Tip	18 G x 1.5 in. BD™ Short Bevel Hypodermic Needle	22 G x 4.0 in. BD™ Short Bevel Hypodermic Needle	22 G x 3.0 in. BD™ Short Bevel Hypodermic Needle	22 G x 2.0 in. BD™ Short Bevel Hypodermic Needle	22 G x 1.5 in. BD™ Short Bevel Hypodermic Needle	Epidural Loss of Resistance (LOR) Syringes	5 mL BD Epilor™ Glass LOR Syringe, Luer Slip	5 mL BD Epilor™ Glass LOR Syringe, Luer Lock, Metal Tip	20 mL BD™ Plastic Syringe, Luer Lock	10 mL BD™ Plastic Syringe, Luer Lock	5 mL BD™ Plastic Syringe, Luer Lock	3 mL BD™ Plastic Syringe, Luer Lock	3 mL BD™ Plastic Syringe, Luer Slip	1 mL BD™ Plastic Syringe, Luer Slip	10 mL BD™ Plastic 3-Ring Control Syringe, Luer Lock	Tetracaine HCL (1%), 2 mL	Lidocaine HCL (1%), 5mL (Skin Wheel)	Lidocaine HCL (1.5%), 20 mL	Sodium Chloride (0.9%), 10 mL (Diluent - LOR)	Drapes	Tray Style (Packaging)	
405901	1T3001	x			x	x	x	x								x	2	x		x		x							P	SD
405904	1T3003		x	x			x											2		x			x					CP	SD	
405905	1T3004		x	x			x											2	x	x					x			P	SD	
405902	2T2235				x		x	x	x	x	x	x	x					x		x			x		x	x	x	P	SD	
405903	2T2385				x		x	x	x	x					x		x				x		x			x		P	SD	
405906	2T0776	x			x		x	x							x			x	x	x					x		x	CP	DD	

Pre-attachments:

3 mL syringe attached to 18 G x 1/2 in. filter needle

Notes:

P is paper drape

CP is clear plastic drape

DD is a BD double-deck tray

SD is a BD single-deck tray with a removable prep tray

For additional information
contact us at 1.888.237.2762
or visit us online at
www.bd.com/anesthesia



1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/anesthesia