



BD™ Spinal Trays

Because your patients trust you to choose the best



Helping all people
live healthy lives

Premium Trays

- Full line of BD™ Whitacre and BD™ Quincke Needle products
- Customized trays available

Precision Components

- High-flow BD™ Whitacre Needles
- Pencil point is designed to spread the dural fibers and help reduce post-dural puncture headaches (PDPH)*
- Precision-formed side hole

BD™ Spinal Trays

Because your patients trust you to choose the best

All Trays Contain:

- 18 G x 1.5 in. BD™ Filter Needle (5 micron), blunt tip
- Prep section including towel, prep cup, prep swabs, and gauze pads
- Povidone iodine topical solution, USP (1% available iodine)
- Fenestrated drape
- Sharps stick pad (not included in BD single-deck configuration)
- Adhesive bandage (not included in BD single-deck configuration)

Material Number	Code Number	26 G x 3.5 in. BD™ Quincke Spinal Needle	25 G x 3.5 in. BD™ Quincke Spinal Needle	22 G x 3.5 in. BD™ Quincke Spinal Needle	25 G x 3.0 in. BD™ Quincke Spinal Needle	25 G x 3.5 in. BD™ Whitacre Spinal Needle, High Flow	22 G x 3.5 in. BD™ Whitacre Spinal Needle	20 G x 1.25 in. Introducer Needle	18 G x 1.25 in. Introducer Needle	27 G x 1.25 in. BD™ Hypodermic Needle	27 G x 0.5 in. BD™ Hypodermic Needle	25 G x 1.5 in. BD™ Hypodermic Needle	25 G x 1.0 in. BD™ Hypodermic Needle	25 G x 1.0 in. BD™ Hypodermic Needle w/ 5Dre-Att: 3mL LL Syringe	22 G x 1.5 in. BD™ Hypodermic Needle	21 G x 1.5 in. BD™ Hypodermic Needle	5 mL BD™ Plastic Syringe, Luer Lock	5 mL BD™ Plastic Syringe, Luer Slip	3 mL BD™ Plastic Syringe, Luer Lock	3 mL BD™ Plastic Syringe, Luer Lock w/ 5Dre-Att: 25 G x 1" Needle	3 mL BD™ Plastic Syringe, Luer Slip	5 mL BD Glaspak™ Syringe, Luer Slip (Spinal)	Tetracaine HCL (1%), 2 mL	Tetracaine HCL Lyophilized, 20 mg	Lidocaine HCL (5%) with Dextrose (7.5%), 2 mL	Bupivacaine (0.75%) with Dextrose (8.25%), 2 mL	Dextrose (10%), 3 mL	Dextrose (10%), 5 mL	Lidocaine HCL (1%), 2 mL (Skin W/Heal)	Lidocaine HCL (1%), 5mL (Skin W/Heal)	Epinephrine (0.1%), 1 mL	Ephedrine Sulfate, 1 mL	Drapes	Tray Style (Packaging)					
405614	1T1122		x										x	x					x	x	x						x			x	x			P	SD				
405616	1T2121		x										x	x					x	x	x						x			x	x			P	DD				
405615	1T1221		x										x	x					x	x		x					x	x			x				P	SD			
405611	1T1321		x										x	x					x	x		x									x					P	SD		
405612	1T2321		x										x	x					x	x		x								x	x					P	DD		
405609	1T2421		x										x	x					x	x										x	x					P	DD		
405618	1T1152	x						x					x	x					x	x		x						x			x	x					P	SD	
405630	1T2151		x					x					x	x					x	x		x						x			x	x					P	DD	
405629	1T1251		x					x					x	x					x	x		x							x			x					P	SD	
405623	1T1351		x					x					x	x					x	x		x									x						P	SD	
405624	1T2351		x					x					x	x					x	x		x									x	x					P	DD	
405621	1T2451		x					x					x	x					x	x											x	x					P	DD	
405635	1T2161	x						x					x	x					x	x		x									x	x					P	DD	
405636	1T2361	x						x					x	x					x	x											x	x					P	DD	
405671	1T2450					x		x					x	x					x	x											x	x					CP	DD	
405607	2T2244			x									x	x						x	x																P	SD	
405607	2T2245			x									x			x					x																P	SD	
405608	2T2246		x					x					x	x						x	x																	P	SD
405610	2T2016			x									x	x						x	x										x							P	SD
405613	2T2218			x									x			x					x	x									x	x						P	SD
405617	2T2119			x									x	x						x	x										x							P	SD
405619	2T2019		x					x					x	x						x	x										x							P	SD
405620	2T2128		x					x					x	x						x	x										x	x	x					P	SD
405622	2T0702		x					x					x							x																		P	DD
405625	2T0146				x			x					x			x				x	x										x							P	SD
405626	2T2018		x					x					x	x						x	x																	P	SD
405627	2T2122		x					x					x	x						x	x																	P	SD
405628	2T0303		x						x				x								x	x																P	SD
405631	2T0251	x						x		x	x		x			x					x	x																P	DD
405632	2T0604	x						x					x			x					x	x																P	DD
405633	2T2125	x						x					x	x						x	x																	P	SD
405634	2T0602	x						x					x			x					x	x																P	DD
405649	2T0045		x				x		x				x				x				x	x																P	SD
405693	2T0086		x		x			x		x			x			x					x	x																P	DD
405694	2T2226		x			x		x					x	x							x	x																P	SD
406701	2T2326		x										x	x							x	x																P	SD

Pre-attachments:

3 mL syringe pre-attached to 18 G filter needle in all codes in BD Single-Deck configurations

Notes:

- P is paper drape
- CP is clear plastic drape
- DD is a BD double-deck tray
- SD is a BD single-deck tray with a removable prep tray
- *Data on file at BD

BD and BD Logo, are trademarks of Becton, Dickinson and Company. © 2007 BD AS1235-0107

For additional information
contact us at 1.888.237.2762
or visit us online at
www.bd.com/anesthesia



1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/anesthesia