

Interlink® System and BD™ Cannulas

Interlink® Locking Cannula Products for use with Interlink System

- 100 per box, 1,000 per case
- 303369** Threaded Lock Cannula
- 303370** Lever Lock Cannula
- 303365** Y Lock Cannula
- 303366** Bare Cannula

Syringes with BD™ Cannula for use with Baxter Interlink®, Abbott LifeShield®, or McGaw SafeLine® Systems

- 100 per box, 800 per case
- 303391** 3ml Syringe with BD Twinpak™ Cannula
- 303346** 3ml Syringe with BD Blunt Plastic Cannula
- 303400** 3ml Syringe with Interlink Cannula
- 303401** 3ml Syringe with Interlink Vial Access Cannula

- 100 per box, 400 per case
- 303347** 5ml Syringe with BD Blunt Plastic Cannula
- 303348** 10ml Syringe with BD Blunt Plastic Cannula
- 303392** 5ml Syringe with BD Twinpak Cannula
- 303393** 10ml Syringe with BD Twinpak Cannula
- 303402** 5ml Syringe with Interlink Cannula
- 303403** 5ml Syringe with Interlink Vial Access Cannula
- 303404** 10ml Syringe with Interlink Cannula

BD Safety-Lok™ Safety Syringes

1ml BD Safety-Lok Insulin Syringe with Permanently Attached Needle

- 100 per box, 500 per case
- 329464** 29g x 1/2 in U-100 BD Ultra-Fine™

1ml BD Safety-Lok Tuberculin Syringe with Permanently Attached Needle

- 100 per box, 500 per case
- 305553** 27g x 1/2 in
- 305554** 25g x 5/8 in

3ml BD Safety-Lok Syringe with Detachable Needle

- 100 per box, 800 per case
- 309592** 25g x 5/8 in **309596** 22g x 1 in
- 309594** 23g x 1 in **309595** 21g x 1-1/2 in
- 309593** 22g x 1-1/2 in

BD Blood Collection Assembly for use with Baxter Interlink®, Abbott LifeShield®, or McGaw SafeLine® Systems

- 25 per box, 200 per case
- 303380** BD Vacutainer® Tube Holder and Luer Adapter with Preattached BD Blunt Plastic Cannula

BD Cannula Products for use with Baxter Interlink®, Abbott LifeShield®, or McGaw SafeLine® Systems

- 100 per box, 1000 per case
- 303345** BD Blunt Plastic Cannula
- 303390** BD Twinpak™ Dual Cannula Device
- 303367** Interlink Vial Access Cannula

BD SafetyGlide™ Syringes and Needles

BD SafetyGlide Hypodermic Needle

- 50 per box, 500 per case
- 305921*** 27g x 5/8 in Needle
- 305901** 25g x 5/8 in Needle
- 305916** 25g x 1 in Needle
- 305902** 23g x 1 in Needle
- 305900** 22g x 1-1/2 in Needle
- 305919*** 22g x 1 in Needle
- 305915** 21g x 1 in Needle
- 305917** 21g x 1-1/2 in Needle
- 305918*** 18g x 1-1/2 in Needle
- * Coming Soon.

BD SafetyGlide Needle with Syringes

- 50 per box, 400 per case
- 1ml BD Syringe with Detachable Needle
- 305903** 25g x 5/8 in
- 305927*** 27g x 5/8 in
- 3ml BD Luer-Lok™ Syringe with Detachable Needle
- 305924*** 25g x 1 in
- 305904** 25g x 5/8 in
- 305905** 23g x 1 in
- 305906** 22g x 1-1/2 in
- 305909** 21g x 1-1/2 in
- 305923*** 22g x 1 in

5ml BD Luer-Lok™ Syringe with Detachable Needle

- 305907** 22g x 1-1/2 in

10ml BD Luer-Lok™ Syringe with Detachable Needle

- 305908** 22g x 1-1/2 in
- * Coming Soon.

BD SafetyGlide Insulin Syringe with Permanently Attached Needle

- 100 per box, 400 per case
- 305930** 1ml Insulin Syringe with 29g x 1/2 in
- 305932** 1/2ml Insulin Syringe with 29g x 1/2 in
- 305935** 3/10ml Insulin Syringe with 29g x 1/2 in
- 305934** 1/2ml Insulin Syringe with 30g x 5/16 in Thin Wall
- 305937** 3/10ml Insulin Syringe with 31g x 5/16 in Thin Wall

BD SafetyGlide Tuberculin Syringe with Permanently Attached Needle

- 100 per box, 400 per case
- 305945** 1ml Tuberculin Syringe with 27g x 1/2 in Regular Bevel
- 305946** 1ml Tuberculin Syringe with 26g x 3/8 in Intradermal Bevel

BD SafetyGlide Allergy Syringe Tray with Permanently Attached Needle

- 25 per tray, 1000 per case
- 305950** 1ml Allergy Tray with 27g x 1/2 in Regular Bevel
- 305951** 1ml Allergy Tray with 26g x 3/8 in Intradermal Bevel

BD Eclipse™ Needles and Syringes

BD Eclipse Needles

- 100 per box, 1,200 per case
- 305767** 25g x 1-1/2 in
- 305761** 25g x 1 in
- 305762** 27g x 1 in
- 305768** 22g x 1 in
- 305763** 22g x 1-1/2 in
- 305765** 21g x 1-1/2 in
- 305764** 21g x 1 in
- 305766** 18g x 1-1/2 in

BD Eclipse Needles with Syringes

- 50 per box, 300 per case
- 1ml BD Luer-Lok™ Syringe with Detachable Needle
- 305778** 30g x 1/2 in
- 305789** 27g x 1 in
- 305780** 25g x 5/8 in

3ml BD Luer-Lok Syringe with Detachable Needle

- 305781** 25g x 5/8 in **305783** 22g x 1-1/2 in
- 305787** 25g x 1 in **305779** 21g x 1 in
- 305782** 23g x 1 in **305784** 21g x 1-1/2 in
- 305788** 21g x 1-1/2 in

5ml BD Luer-Lok Syringe with Detachable Needle

- 305785** 22g x 1 in

10ml BD Luer-Lok Syringe with Detachable Needle

- 305786** 22g x 1-1/2 in

BD Integra™ Needles and Syringes

BD Integra 3ml Syringe with Detachable Needle

- 100 per box, 400 per case
- 305269** 3ml Syringe with 25g x 5/8 in Needle
- 305270** 3ml Syringe with 25g x 1 in Needle
- 305271** 3ml Syringe with 23g x 1 in Needle
- 305272** 3ml Syringe with 22g x 1-1/2 in Needle
- 305273** 3ml Syringe with 21g x 1 in Needle
- 305274** 3ml Syringe with 21g x 1-1/2 in Needle

BD Integra Needles

- 100 per box, 1,000 per case
- 305310** 25g x 5/8 in Retracting Needle
- 305311** 25g x 1 in Retracting Needle
- 305312** 23g x 1 in Retracting Needle
- 305313** 22g x 1-1/2 in Retracting Needle
- 305275^{1,2}** 18g x 1-1/2 in Blunt Fill Needle
- 305833^{1,2}** 18g x 1-1/2 in Blunt Filter Needle

- 1 – 18g Blunt Fill and Blunt Filter Needles have a non-retracting needle
- 2 – Packaged 50 Box/500 Case

BD Integra 1ml Syringe with Permanently Attached Needle

- 100 per box, 800 per case
- 305282** 1ml Insulin Syringe with 29g x 1/2 in
- 305298** 1ml TB Syringe with 27g x 1/2 in
- 305279** 1ml TB Syringe with 26g x 3/8 in
- 305278** 1ml Syringe with 25g x 1 in
- 305277** 1ml Syringe with 25g x 5/8 in

BD Integra 3ml Syringe Only

- 100 per box, 400 per case
- 305283** 3ml Syringe Only

BD PosiFlush™ Prefilled Syringes

BD™ Prefilled Flush Syringes Are Available In A Wide Variety of Sizes and Solutions To Meet All Your Flush-Protocol Requirements.

Product #	Description	Units/Box	Boxes/Case	Units/Case
BD PosiFlush™ Syringes Filled With Saline				
306570	BD saline XS 3mL	30	8	240
306571	BD Saline XS 5mL	30	8	240
306572	BD Saline XS 10mL	30	8	240
306554	BD Saline SP 3mL	30	8	240
306555	BD Saline SP 5mL	30	8	240
306556	BD Saline SP 10mL	30	8	240

New! BD PosiFlush Heparin Lock Flush Syringes coming soon.



Hypodermic Product Catalog

The Basis for Selection of Your Medication Delivery Device Needs



BD Medical
2771 Bristol Circle
Oakville, ON L6H 6R5
www.bd.com

All products are Latex-Free, Preservative-Free

NOTE: Bulk non-sterile syringes in bulk packaging are also available.

For more information please call: 1.800.268.5430

Interlink is a registered trademark of Baxter International, Inc.
Lifeschild is a registered trademark of Abbott Laboratories.
SafeLine is a registered trademark of B. Braun, Inc.
BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2005 BD.
Printed in Canada BDM0094- CE 03/05

BD™ Insulin Syringes

3/10ml BD Lo-Dose™ Insulin Syringe with Permanently Attached Needle

Self-Contained, 10 ea. per bag, 200 per box, 1,600 per case
309300 28g x 1/2 in U-100 BD Micro-Fine™ IV
309301 29g x 1/2 in U-100 BD Ultra-Fine



1/2ml BD Lo-Dose Insulin Syringe with Permanently Attached Needle

Blister Package. 100 per box, 500 per case
329461 28g x 1/2 in U-100 BD Micro-Fine IV (Orange)
 200 per box, 1,600 per case
309305 28g x 1/2 in U-100 BD Micro-Fine IV (Orange)
309306 29g x 1/2 in U-100 BD Ultra-Fine



NOTE: All of the syringes in this section are in an Easy-Peel Blister Package unless otherwise noted as Self-Contained Package.

BD™ Allergy Syringes

BD SafetyGlide™ Allergy Syringe Tray with Permanently Attached Needle

25 per tray, 1,000 per case
305950 1ml Allergy Tray with 27g x 1/2 in Regular Bevel
305951 1ml Allergy Tray with 26g x 3/8 in Intradermal Bevel



1ml BD Allergist Tray with Permanently Attached Needle

25 syringes per tray, 1,000 per case
305537 26g x 1/2 in Regular Bevel
305538 26g x 1/2 in Intradermal Bevel
305539 26g x 3/8 in Intradermal Bevel
305540 27g x 1/2 in Regular Bevel
305541 27g x 3/8 in Intradermal Bevel
305542 27g x 3/8 in Regular Bevel



1ml BD Insulin Syringe with Permanently Attached Needle

Self-Contained Package. 200 per box, 1,600 per case
309309 28g x 1/2 in U-100 BD Micro-Fine IV (Orange)
309311 29g x 1/2 in U-100 BD Ultra-Fine
 Blister Package. 100 per box, 500 per case
329424 28g x 1/2 in U-100 BD Micro-Fine IV (Orange)
329420 28g x 1/2 in U-100 BD Micro-Fine IV Single Unit Graduation



1ml BD Insulin Syringe with Detachable Needle

Blister Package. 100 per box, 800 per case
329651 25g x 5/8 in U-100 (Orange/Blue)
329652 26g x 1/2 in U-100 (Orange/Tan)
329622 25g x 1 in U-100 (Orange/Blue)



1ml BD Insulin Syringe Only


Blister Package. 100 per box, 800 per case
329650 U-100 Slip Tip (Orange)



BD™ Tuberculin Syringes

1/2ml BD Tuberculin Syringe with Permanently Attached Needle

100 per box, 500 per case 1/100ml
305620 27g x 1/2 in Regular Bevel



1ml BD Tuberculin Syringe Only

100 per box, 800 per case 1/100ml
309602 Slip Tip



BD™ Syringes & Needles (sterile, single-use)


1ml BD Syringe/Needle Combination, Slip Tip

100 per box, 800 per case 1/100ml
309597 26g x 5/8 in Sub-Q



1ml BD Luer-Lok™ Syringe/Needle Combination

100 per box, 800 per case 1/100ml
309637 20g x 1 in



1ml BD Luer-Lok™ Syringe Only

100 per box, 800 per case 1/100ml
309628 BD Luer-Lok™



1ml BD Tuberculin Syringe with Detachable Needle, Slip Tip

100 per box, 800 per case 1/100ml
309623 27g x 1/2 in
309625 26g x 3/8 in Intradermal Bevel
309626 25g x 5/8 in
309624 21g x 1 in



3ml BD Syringe/Needle Combination, BD Luer-Lok™ Tip

100 per box, 800 per case 1/10ml
309582 25g x 1-1/2 in
309581 25g x 1 in
309570 25g x 5/8 in
309571 23g x 1 in
309574 22g x 1-1/2 in
309573 22g x 1-1/4 in
309572 22g x 1 in
309569 22g x 3/4 in
309577 21g x 1-1/2 in
309575 21g x 1 in
309579 20g x 1-1/2 in
309578 20g x 1 in
309580 18g x 1-1/2 in
305060 18g x 1-1/2 in Blunt Fill Needle



3ml BD Syringe Only

100 per box, 800 per case 1/10ml
309585 BD Luer-Lok™ Tip
309586 Slip Tip



5ml BD Syringe/Needle Combination, BD Luer-Lok™ Tip

100 per box, 400 per case 1/5ml
309630 22g x 1 in
309631 22g x 1-1/2 in
309632 21g x 1 in
309633 21g x 1-1/2 in
309634 20g x 1 in
309635 20g x 1-1/2 in



5ml BD Syringe Only

100 per box, 400 per case 1/5ml
309603 BD Luer-Lok™ Tip
301603 Slip Tip



BD™ Syringes & Needles (sterile, single use)

10ml BD Syringe/Needle Combination, BD Luer-Lok™ Tip

100 per box, 400 per case 1/5ml
309640 22g x 1-1/2 in
309642 21g x 1 in
309643 21g x 1-1/2 in
309644 20g x 1 in
309645 20g x 1-1/2 in



10ml BD™ Syringe Only

100 per box, 400 per case 1/5ml
309604 BD Luer-Lok™ Tip
301604 Slip Tip
305482 Eccentric Tip



BD Pharmacy Convenience Paks

Tray Contents:
309701 1ml Luer Slip Tip 25 syringes only per sterile tray, 12 trays per case of 300
309702 3ml BD Luer-Lok Tip 25 syringes only per sterile tray, 12 trays per case of 300
309703 5ml BD Luer-Lok Tip 25 syringes only per sterile tray, 12 trays per case of 300
309605 10ml BD Luer-Lok Tip 20 syringes only per sterile tray, 12 trays per case of 240
305617 20ml BD Luer-Lok Tip 10 each Syringes only per sterile tray, 12 trays per case of 120
305618 30ml BD Luer-Lok Tip 10 each Syringes only per sterile tray, 12 trays per case of 120
309680 60ml BD Luer-Lok Tip 20 each Syringes only per sterile tray, 6 trays per case of 120



BD Luer-Lok Sterile Syringe Tip Cap

200 per box, 2000 per case
305819 BD Syringe Tip Cap




20ml BD Syringe Only

40 per box, 160 per case 1ml, 3/4 oz. in 1/8 oz.
309661 BD Luer-Lok Tip
301625 Slip Tip
 120 per box, 480 per case
300613 Eccentric Tip



30ml BD Syringe Only

40 per box, 160 per case 1ml, 1 oz. in 1/4 oz.
309650 BD Luer-Lok Tip
309651 Slip Tip



60ml BD Syringe Only

40 per box, 160 per case 1ml, 2 oz. in 1/2 oz.
309654 Slip Tip
 40 per box, 160 per case
309653 BD Luer-Lok Tip
 60 per box, 240 per case
300866 Eccentric Tip



2oz BD Catheter Tip Syringe Only (with Tip Shield)

40 per box, 160 per case 1ml, 2 oz. in 1/2 oz.
309620 Catheter Tip



BD PrecisionGlide™ Needles (sterile, single-use)

BD™ Regular Bevel Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305106	30g x 1/2 in	
305128	30g x 1 in	
301629	27g x 1-1/2 in	
305109	27g x 1/2 in	
305136	27g x 1-1/4 in	
305110	26g x 3/8 in	Intradermal Bevel
305111	26g x 1/2 in	
305122	25g x 5/8 in	
305124	25g x 7/8 in	
305125	25g x 1 in	
305127	25g x 1-1/2 in	
305143	23g x 3/4 in	
305145	23g x 1 in	
305120	23g x 1-1/4 in	
305155	22g x 1 in	
305156	22g x 1-1/2 in	
305165	21g x 1 in	
305167	21g x 1-1/2 in	
305129	21g x 2 in	
305175	20g x 1 in	
305176	20g x 1-1/2 in	
305186	19g x 1 in	Thin Wall
305187	19g x 1-1/2 in	Thin Wall
305195	18g x 1 in	
305196	18g x 1-1/2 in	
305197	16g x 1 in	
305198	16g x 1-1/2 in	
305115	26g x 5/8 in	Sub-Q
305193	23g x 1 in	Thin Wall IM
305194	23g x 1-1/2 in	Thin Wall IM
305185	18g x 1-1/2 in	Thin Wall Fill
305177	21g x 1 in	Thin Wall IV

BD™ Blunt Fill Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305180	18g x 1-1/2 in	
305211	18g x 1-1/2 in	Filter (5 Micron)

BD™ Short Bevel Needles

100 per box, 1,000 per case

Reorder Number	Needle Size	Description
305159	22g x 1-1/2 in	
305178	20g x 1 in	
305179	20g x 1-1/2 in	
305188	19g x 1 in	Thin Wall
305189	19g x 1-1/2 in	Thin Wall
305199	18g x 1-1/2 in	

Types of needle walls & bevels

The Thin Wall needle is thinner than a Regular Wall needle. The inner diameter is bigger without any increase in the outside diameter. For example, a 21g Thin Wall needle has the same outside diameter as a 21g Regular Wall, but has the same inner diameter as a 20g needle. Thin Wall allows the same flow volume as a Regular Wall needle one gauge larger. Thin Wall is available in 18, 19, 21, 22, and 23 gauge sizes.


Regular Wall



Thin Wall



Regular Bevel



Short Bevel



Intradermal Bevel

BD™ Sharps Disposal Collectors

Point-of-Use Product

Size (L)	Package Quantity	#/Case	Description
300475	5.1	12	Horizontal Entry
300974	5.1	20	Horizontal Entry
303050	7.6	10	Horizontal Entry
303051	11.3	10	Horizontal Entry

Nestable Product

Size (L)	Package Quantity	#/Case	Description
300450	3.1	36	Inner Liner (yellow)
305469	3.1	36	Inner Liner (clear)
300439	7.6	24	Funnel Cap
300443	13.2	20	Open Cap
300434	2.7	12	Point-first Vertical Entry
300483	2.7	12	X-Large Funnel Cap
300436	n/a	1	Physician Office Safety Kit
300438	n/a	6	Physician Kit Refills
300988	13.2	20	Funnel Cap
305555	7.6	24	Open Cap

Chemotherapy Product

Size (L)	Package Quantity	#/Case	Description
300459	6.6	12	Plug cap
300162	10.3	12	X-large Plug Cap
300161	18.7	8	Plug Cap

Phlebotomy Product

Size (L)	Package Quantity	#/Case	Description
367220	13.2	20	Universal Unwinder
367224	3.0	30	Universal Unwinder
367216	1.0	60	Universal Unwinder

One Piece Product

Size (L)	Package Quantity	#/Case	Description
300460	1.4	36	Funnel Cap
300462	1.4	36	Plug Cap
300466	3.1	24	Funnel Cap
300452	10.3	12	X-large Plug Cap

Extra Large Product

Size (L)	Package Quantity	#/Case	Description
305579	34	8	XL Sharps Collector
305647	72	8	XL Sharps Collector

Accessories Accessory

Description	Package Quantity
305447 Locking Wall Bracket	10
305017 Patient Room Cabinet	1
305409 Bracket	10
305494 Security Padlock	12
305475 Wall Cabinet	12
305485 Bracket for Nestable Containers	12
305497 Quick Release Bracket	12
305495 X-large Security Bracket	12
305451 Tray Bracket	12
305647 9 gal. Foot Operated Trolley	1
305442 Base Stand Floor Bracket	5



30G

27G

26G

25G

24G

23G

22G

21G

20G

19G

18G

16G

BD Conventional Needle & Syringe Color Chart