

BD infusion solutions for TPN

Adapting treatments, not lifestyle.

Malnutrition is a major clinical condition that can be caused by, among other things, gastrointestinal dysfunctions, surgery or cancer. Patients with wide ranging medical conditions rely on the use of total parenteral nutrition (TPN). However, a TPN pump can compromise the patient's quality of life by giving them less independence and mobility.



The BD BodyGuard[™] is a small, intuitive and lightweight multi-therapy ambulatory infusion pump that helps patients achieve their desired independence.

Intuitive and easy to program, helping reduce patient/ nurse training:

- Large color screen and onscreen guided messages
- Easily programmable in a few steps
- Delayed start functionality
- Last program stored functionality
- Priming via the pump

Fully fitted for TPN:

- Programmable mode options that include ramp-up/ ramp-down and continuous infusion modes
- Wide range of accessories specific for TPN adult, pediatric use and homecare use
- Range of sets specifically created for TPN administration with a 1.2 micron filter
- Silent operation


Advanced safety features to help ensure patient safety

- Lock mode (screen lock, program lock, rate lock)
- 3 level access codes to help prevent unauthorised program changes
- Extensive Drug Library and fully configurable drug limits
- Audible alarms including occlusion, upstream/ downstream and end of infusion
- Sets with built-in filters and a manual priming valve
- Free flow and back flow protection
- Onscreen pressure bar graph





What you need to start your TPN infusion



Contact your BD representative for a demonstration.

2L BD BodyGuard[™] Bag

1. ESPEN guidelines on chronic intestinal failure in adults, January 2016

BD Switzerland Sàrl, Terre Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland.



bd.com

BD, the BD logo and BodyGuard are trademarks of Becton, Dickinson and Company or its affiliates. ©2020 BD. All Rights Reserved. 0000CF03484 Issue 2.