



## Alaris® Enteral syringe pump

Designed to deliver nutrition efficiently, the Alaris® Enteral syringe pump includes a large quick-to-read screen, intuitive keyboard, alarm beacon and familiar user interface. The Alaris® Enteral syringe pump recognises a wide range of enteral syringes and is clearly identified with its purple labelling for fast identification of enteral use.

# Alaris® Enteral syringe pump



Designed to deliver nutrition efficiently, the Alaris® Enteral syringe pump includes a large quick-to-read screen, intuitive keyboard, alarm beacon and familiar user interface. The Alaris® Enteral syringe pump recognises a wide range of enteral syringes and is clearly identified with its purple labelling for fast identification of enteral use.



## Alaris® Enteral syringe pump Specifications

### Rate Range

5ml syringe	–	0.1 to 150ml/h
10ml syringe	–	0.1 to 200ml/h
20ml syringe	–	0.1 to 200ml/h
30ml syringe	–	0.1 to 200ml/h
50ml syringe	–	0.1 to 200ml/h

### Volume Delivered

0.0ml - 9990ml

### Volume To Be Delivered range

0.1ml - 1000ml, 1 min - 24hrs

### Bolus Rate

5ml syringe	–	10 to 150ml/h
10ml syringe	–	10 to 200ml/h
20ml syringe	–	10 to 200ml/h
30ml syringe	–	10 to 200ml/h
50ml syringe	–	10 to 200ml/h

### Max Purge Rate

500ml/h (limited to size of syringe)

### Max Purge Volume Range

Purge volume range is 2ml

### Volume to be delivered over time

Pre-defined volume can be set to deliver over pre-defined time.

### Battery Type

Rechargeable sealed NiMH. Automatically charges when the pump is connected to AC power.

### Battery Life

Mean Time To Battery Empty from fully charged at 5ml/h & 23°C under normal conditions is 6 hours.  
(95% lower confidence interval of 5 hours 50 minutes)

### Recharge Time

Charging takes 2½ hours from discharge to 90% charge.

### Pumping Pressure Limits

11 Levels L0-L10 (50mmHg - 1000mmHg)

### Alarms, Warnings and Prompts

Audible alarm and display prompts with comprehensive protection such as:

- Drive disengaged
- Check syringe
- Near end of delivery (NEOD)
- VTBD Done
- Internal malfunction
- Occlusion
- Battery low
- End of Delivery (EOD)
- Attention (nurse callback)

### AC Power Requirements

115-230 VAC, 50/60Hz. 20VA (nominal)

### Dimensions

310mm(w) x 121mm(h) x 200mm(d)

### Weight

2.7 Kg

	5ml	10ml	20ml	30ml	50ml	60ml
CareFusion Enteral	0000ME00870	0000ME00871	0000ME00872	0000ME00873	✗	0000ME00874
Pentaferite Enteral	2022590	2022690	2022790	✗	✗	2022990
Vygon Enteral	✗	1015.102	1015.212	✗	✗	1015.602
Vygon A-VY Enteral	✗	✗	1015.213	✗	✗	1015.603
Medicina B.Tip	✗	✗	✗	✗	✗	Bladder Tip PE60B
Medicina Enteral	PE05	PE10	PE20	✗	✗	PE60
Nutricair	✗	SE10	SE20	✗	SE50	✗
Terumo Enteral	✗	✗	✗	SS+30EO	SS+50EO/C	✗

CareFusion International Headquarters  
A-One Business Centre  
Zone d'activités Vers-la-Pièce no 10  
1180 Rolle  
Switzerland  
+41 21 556 30 00 tel  
+41 21 556 30 99 fax

CareFusion South Africa  
Unit 2 Oude Molen Business Park  
Oude Molen Road  
Ndabeni  
Cape Town 7405  
South Africa  
+27 21 510 7562 tel  
+27 21 510 7567 fax

CareFusion UK Ltd  
The Crescent, Jays Close, Basingstoke  
Hampshire RG22 4BS  
United Kingdom  
0800 917 8776 freephone  
+44 (0)1256 330860 fax

CareFusion  
PO Box 5527  
Dubai  
United Arab Emirates  
Tel: +9714 2822 842  
Fax: +9714 2822 914

© 2010 CareFusion Corporation or one of its subsidiaries. All rights reserved.  
Alaris is a registered trademark of CareFusion Corporation or one of its subsidiaries. 0000CF00567 Issue 2