



Oncology Ambulatory Infusions

Efficient. Flexible. Safe.





Meeting your medical oncology treatment needs with our ambulatory infusion solutions.

At BD, we understand that oncology is a busy and time-sensitive practice. Ensuring low patient waiting times and higher satisfaction rates are imperative. As the drugs used in oncology treatment are expensive and often hazardous, drug wastage reduction and healthcare professional safety is crucial.

From preparation to administration, our holistic and efficient ambulatory infusion solution help ensure that:

- Patients and caregivers stay safe and confident during therapy
- Drug wastage is reduced, and workflow efficiency is increased



Increased efficiency and improved patient experience.

Efficient drug and waste management is key, along with delivering fast patient treatment times. The BD BodyGuard™ Duo dual channel infusion pump and BodyGuard™ MicroSets help ensure that time, cost, and space savings are made without sacrificing patient experience or safety, so healthcare professionals can keep up with whatever each day brings.

Improved efficiency allows more patients to be served, whilst reducing cost per treatment and drug wastage. The compact size and intuitive setup of our pump means that practitioners can enable treatment as quickly and effectively as possible, leaving them more time to focus on delivering excellent patient care. The dual-channel pump and microbore set offering:

- Aids patient mobility
- Enables concurrent flow of medications (two different infusions delivered at the same time): meaning better adherence to medication therapy and fewer pumps required
- Reduces flush times between chemotherapy treatments when compared to industry standard diameter tubing



Greater flexibility ensures greater control.

Oncology is fast-paced, and best practice continues to evolve and shift. The BodyGuard™ Duo offers the features and capabilities you need to respond rapidly - now and into the future. An ambulatory dual channel pump, alongside our flexible service and training offer, will help you to stay nimble and in charge as the oncology landscape changes.

In the fast-paced environment of oncology daycare units, which are inherently susceptible to delays, greater flexibility assists you in taking back control.

The intuitive pump, one of the smallest smallest and lightest available, delivers a variety of therapies in a single device (so no need for additional pumps). The selection of low priming

volume administration sets, that are multi-channel and multi-therapy, enable you to address almost all therapeutic requirements while reducing time and drug wastage.

The mono and adjacent therapy solutions in BD's comprehensive portfolio of IV sets for oncology (BD BodyGuard™ Microsets) include primary and secondary systems; low sorbing materials; SmartSite™ connectors to help reduce accidental exposure hazards; and low wastage, helping to support you in the delivery of efficient, safe, cost-effective treatment for more patients.

We also provide flexible service programs and training which make it easier for you to implement our solutions and stay on top of BD technology.

Be confident in performance and safety.


BD has been leading the way in patient safety for over 100 years. The secure delivery of hazardous drugs in an oncology setting, and improving the well-being of patients and practitioners, continues to be our priority.

The BD BodyGuard™ MicroSets closed system help ensure safety during treatment preparation and administration.

BD's broad portfolio of infusion sets helps to deliver free flow protection, reduce the risk of exposure to dangerous drugs, and improve workflow efficiency.

Our training and servicing packages provide additional confidence in pump performance and safety – treatment after treatment. Ensuring a holistic safety approach to oncology treatment.





A total care solution and end-to-end support – all day, every day.

Oncology treatment is lifesaving, but it's also complicated and toxic. You're only as strong as your weakest link. That's why BD supports you in everything from prescription to administration, all the way through to training, service plans and clinical consulting. We work with you to help you get it right, day after day, year after year, whatever the care setting.

Our unified approach to training and support helps you to reduce the complexity of your treatment offering. Empowering your practitioners to follow key best practices simply and effectively, offers peace of mind, and a potential reduction in total treatment delivery costs.

We can provide service and upkeep contracts, as well as complementary products and services such as the BD PhaSeal™ closed system transfer devices, BD Cato™ IV preparation, training, service plans, and clinical consulting.


In summary, we can give you all the support you need to deliver oncology treatment efficiently, flexibly, safely, and cost-effectively.

From prescription to administration

To help you achieve greater efficiency, flexibility, and safety, we offer innovative solutions in every step of your medical oncology treatment.

The BD Cato™ IV Prep system is a guided, gravimetrics-based IV workflow management system designed to help support your pharmacy's compounding operations for all preparation types, including IVs, hazardous and non-hazardous preps, sterile and non-sterile preps and oral syringes.

The BD PhaSeal™ leading Closed System Transfer Device (CSTD) has been demonstrated to help prevent exposure to hazardous drugs, from drug preparation to IV administration.




BD BodyGuard™ IV MicroSets


Alongside our leading medication preparation devices in the oncology environment, we offer a range of complementary BD BodyGuard™ IV sets that provide the following benefits:

- Multi-way and secondary extension sets help easily and safely mix hazardous medications in the pharmacy with built in multi needle-less SmartSite™
- Administer drugs concurrently while helping limit drug wastage
- Low sorbing, PE-lined tubing helps minimize the interaction of tube material with hazardous medications
- 1.6mm internal diameter Microbore tubing, helps reduce:
 - internal volume and resulting medication waste
 - mid therapy medication flush times
- With or without built-in 0.2µ filtration
- Ability to deliver blood without risk of red blood cell damage²
- Anti-syphon valves as a standard provide independent free-flow and reflux protection
- Light resistant version for UV sensitive medications
- Drop sensor compatible to help prevent introduction of air into the IV set


Multi-way and Secondary IV Sets



Primary IV Sets



True Dual Channel



With the BD BodyGuard™ Duo pole mount charger, you have the benefit of a far view display.

To help you stay in touch with patients' infusion status, the most valuable infusion parameters toggle and can be viewed from distance as well as wider angle.



To order or learn more about how BD can assist you in your efforts to bring this valuable solution to more patients, contact your BD sales representative.

1. Closed systems are achieved using specific BD products. Contact your local BD representative to learn more.
2. Restriction may apply; refer to pump Direction For Use

BD Switzerland Sàrl, Terre Bonne Park – A4 Route de Crassier 17, 1262 Eysins, Switzerland

[bd.com](https://www.bd.com)

BD, the BD logo, BodyGuard, SmartSite, PhaSeal, Cato, are trademarks of Becton, Dickinson and Company or its affiliates.
©2020 BD All Rights Reserved. 0000CF04871 Issue 1

