

UK Medical

It's the attention in detail that makes the difference!

For over 20 years, UK Medical has been supplying medical products to the NHS and private healthcare sector. We select products that add value and influence patient care pathways for the better.

In an industry that demands high standards, good quality products and a professional service, UK Medical strives to exceed expectations. We believe we are different. We hope that shines through.

ELLA-CS

Founded in 1991, ELLA-CS specialises in the design and manufacture of gastrointestinal stents. Working in close partnership with leading global clinicians and medical institutions the company develops innovative products to meet the challenges of modern gastrointestinal stenting.

A central European location and flexible approach enables ELLA to respond quickly to clinician's needs, something that the company prides itself on. ELLA is responsible for the design and manufacture of some truly ground-breaking technologies including the world's first commercially available biodegradable GI stents, the innovative Danis stenting system and unique ELLA Stent Extractor.

ELLA-CS, Czech Republic


www.ella-cs.eu



UKMED!CAL

stenting@ukmedical.com
www.ukmedical.com

UK Medical Limited
Albreda House,
Lydgate Lane, Sheffield,
S10 5FH United Kingdom
Tel +44 (0)114 268 8880
Fax +44 (0)114 268 8881



Gastrointestinal Stents

www.ukmedical.com



S-Ella-Bro-v1-0911
Date of preparation: Sept 2011

Gastrointestinal Stents

Our comprehensive range of ELLA gastrointestinal stents combines market leading technology and innovative product design to meet the challenges of modern gastrointestinal stenting.

ELLA stents are available in a wide range of diameters, lengths and delivery systems, providing flexible therapeutic solutions to a broad spectrum of digestive disorders.


Nitinella Plus

Nitinella Plus is a highly flexible nitinol biliary SEMS with reduced shortening andatraumatic ends. Available either uncovered or with full or partial silicone covering and optional transhepatic or endoscopic delivery system.





Enterella Pyloroduodenal

Enterella Pyloroduodenal is a highly flexible uncovered nitinol SEMS with excellent radial force, optional TTS or OTW delivery system and proximal repositioning loop.



Enterella Colonic

Enterella Colonic is a highly flexible nitinol SEMS available covered or uncovered, with optional OTW or TTS delivery system proximal and repositioning loop.




Flexella Plus

A proximal release covered nitinol oesophageal SEMS designed for accurate stent placement in the upper oesophagus.



HV Plus

HV Plus is a removable nitinol oesophageal SEMS with full silicone covering, unique anti-migration design and optional polyethylene anti-reflux valve.




BD Biodegradable

BD Biodegradable is the first fully biodegradable oesophageal stent available on the market and is indicated as an alternative to conventional treatment of benign oesophageal strictures (including peptic, anastomotic and caustic) and achalasia refractory to standard therapy.




Boubella

Boubella is a fully covered stainless steel oesophageal SEMS with segmented body design for excellent radial force, proximal flared segment to reduce the risk of migration and optional anti-reflux valve.



Danis

Danis is a fully silicone covered nitinol SEMS uniquely designed to tamponade bleeding oesophageal varices through direct compression.



Ella Stent Extractor

A unique system designed for safeatraumatic removal of ELLA HV Plus, Flexella Plus and Danis stents

