


SmartSite™

Efficient. Flexible.
Standardised.

A widely used technology offering a leak-free system for chemotherapy when used in combination with Texium® devices.

The SmartSite™ connector offers a vial to vein solution, from drug preparation to administration, whilst also enabling connector standardisation with BD infusion pumps.



- 1 Split septum.**
When needleless systems are used, a split septum may be preferred.¹
- 2 Easy to disinfect.**
A smooth easy-to-scrub surface helps reduce the risk of bacterial ingress.²

- 3 Easy to flush.**
The straight fluid pathway provides laminar flow and easy flushing.²
- 4 Compact, lightweight design.**
Suitable for neonate and paediatric use. Easy handling for healthcare workers.

1. O'Grady N, Alexander M, Burns L et al. Guidelines for the Prevention of Intravascular Catheter- Related Infections. Clin Infect Dis. 2011 May 1; 52(9): 197–200.
2. Jarvis W. Choosing the Best Design for Intravenous Needleless Connectors to Prevent Healthcare-Associated Bloodstream Infections. Infection Control Today. August 2010, Vol. 14, No. 8.

BD Switzerland Särl, Terre Bonne Park – A4 Route de Crassier 17, 1262 Eysins, Switzerland

bd.com