

S 1130.1
Revised: September 2002

 **BD Difco™ Trypsin**

INTENDED USE

Difco Trypsin is a stable, filter-sterilized, standardized and desiccated tryptic enzyme used in a variety of procedures, including tissue culture, requiring proteolytic enzymatic activity.

When rehydrated according to directions, **Difco** Trypsin has a tryptic activity equivalent to a 5% solution of Trypsin 1:250, Difco.

Precautions: For Laboratory Use

WARNING: HARMFUL BY INHALATION AND IF SWALLOWED. IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN. May cause sensitization by inhalation and skin contact. Do not get in eyes, on skin, on clothing. Wash thoroughly after handling.

FIRST AID: If swallowed, rinse mouth and throat with water. Call a physician. In case of contact, immediately flush eyes with plenty of water for at least 15 min. Wash exposed skin with soap and water. If irritation persists, seek medical attention.

Storage: Store **Difco** Trypsin at 2 – 8°C.

Rehydration: Rehydrate with 10 mL sterile purified water.

AVAILABILITY

Cat. No. Description

215310	Difco Trypsin, 6 x 10 mL.
215320	Difco Trypsin, 12 x 10 mL.

BD and BD Logo are trademarks of Becton, Dickinson and Company. Difco is a trademark of Difco Laboratories, a subsidiary of Becton, Dickinson and Company. © 2002 BD.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA
800-638-8663