

MARSTERS Incubator #425380

#	<u>DESCRIPTION</u>	<u>CATALOG NUMBER</u>	<u>REFERENCE NUMBER</u>
1	Switch assembly	421814	42538002
2	Lamp Assembly	421817	42538005
3	Right angle thermometer	421818	42538006
4	Heater assembly	421832	42538021
5	Feet (4/pk)	421825	42538013
6	Inserts (5/pk)	421826	42538014
9	PC Board assembly	421829	42538017
10	Thermistor PC board	421830	42538018
11	Potentiometer assembly	421831	42538019

Discontinued parts

8	Top plate	discontinued 8/02 421828	42538016
N/A	Pilot Light	discontinued 8/02 421820	42538008
N/A	Bulb (therm light)	discontinued 8/02 421819	42538007
7	Plugs (5/pk)	discontinued 8/02 421827	42583015




Figure 4. MARSTERS Incubator Illustrated Spare Parts Breakdown.