

Accuracy, Sensitivity, Speed and Convenience

The BD *CMVscan*[™] Latex Agglutination Test accurately detects both IgG and IgM antibodies in serum or plasma.



Only the BD *CMVscan* test gives you clear antigen agglutination readings directly from serum or plasma.* The rapid, simple procedure remains a laboratory favorite. And the convenient kit suits your laboratory work flow perfectly.

Qualitative results in only 4 main steps. 10-minute CMV antibody detection confirms exposure to cytomegalovirus.

Quantitative results in just 5 main steps. 15-minute CMV antibody titer results accurately track infection.

* Li, Y.-C., and I. Rutherford. Comparison of Two Serological Methods for the Detection of Antibody to Cytomegalovirus. Cleveland Clinic Foundation. Abstr. Annu. Meet. Am. Soc. Microbiol. 1987, C-46, p. 331.

Feature	Benefit
Detects both IgG and IgM; detects low levels of antibody.	Early detection; increased sensitivity.
Sensitivity 99%; specificity 97%.	Confidence in results.
Results in 10-15 minutes.	STAT testing capabilities; free up techs for other duties.
Latex agglutination procedure.	Familiar with technology; easy transition.
Qualitative and quantitative tests use same test components.	Less inventory to manage; fewer variables; fewer storage requirements.
No serum pretreatment or antigen titration required.	Saves time; fewer variables; eliminates effects of non-specific factors.
Ready to use reagents.	Saves time; ensures consistency
No expensive instrumentation.	No money spent on capital equipment; no service contracts; not dependent on instrumentation.

Ordering Information

Product Name	Product/Description	Catalog No.
<i>CMVscan</i> Test Kit	30-determination kit complete with all reagents, controls, buffer, supplies and 4 test cards.	255126
<i>CMVscan</i> Test Kit	100-determination kit complete with all reagents, controls, buffer, supplies and 5 test cards.	255001
<i>CMVscan</i> Test Kit	500-determination kit complete with all reagents, controls, buffer, supplies and 20 test cards.	255205



Indispensable to human health

BD Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152-0999 USA
800.638.8663
www.bd.com/clinical