

A Powerful Team — BACTEC™ MGIT™ 960 & EpiCenter™™



Indispensable to
human health

**Multiplying
Your Options
in Mycobacteria
Data Management**


Phoenix™ 100⁺
BACTEC™ 9000 Series
BDProbeTec™ ET



BACTEC™ MGIT™ 960 & EpiCenter™


TB Method Comparison Capability

- Calculates average time-to-detection for individual methods, assuring that CDC guidelines¹ are met for recovery times
- Now you can compare automated methods with classical methods
- Enables bi-directional interface to an LIS




Recovery Rate Reporting

- Store data for 5 years or longer
- Conduct longitudinal tracking of individual TB patients
- Search for shifts in organism and susceptibility patterns



Enables Long Term Data Storage with Historical Trending Capability

- Create customizable reports
- Monitor trendshifts of different organism populations
- Generate monthly trend report



¹Tenover, F., Crawford, J., Huebner, R., Geiter, L., Horsburgh, C., Good, R. 1993. Guest Commentary. The Resurgence of Tuberculosis: Is Your Laboratory Ready? J. Clin. Microbiol. 31:767-770.



7 Loveton Circle
Sparks, MD 21152-9956 USA
800.638.8663
www.bd.com/biosciences/diagnostic

11 rue Aristide Bergès
38800 Le Pont de Claix, France
Tel: 33.4.7668.3636

2771 Bristol Circle
Oakville, Ontario
Canada L6H 6R5
Tel: 800.268.5430

30 Tuas Avenue 2
Singapore 639461
Tel: 65.8610.633

Monte Pelvoux 111, 9th Floor
Col. Lomas de Chapultepec
11000 México D.F.
Tel: 52.5.237.1200

Rua Alexandre Dumas 1976
04717-004 São Paulo, S.P. Brazil
Tel: 55.11.5185.9833

Akasaka DS Building
5-26 Akasaka 8-chome
Minato-ku
Tokyo, 107 Japan
Tel: 81.24.593.5405

BD, Logo and all trademarks are the property of Becton, Dickinson and Company. ©2001. 1-2423 January 2001. Printed in USA.