



Influenza A?

- Positive?
- Rapid Result?
- Negative?
- Which Treatment?

Influenza B?

- Positive?
- Rapid Result?
- Negative?
- Which Treatment?



Indispensable to
human health

One Test, Two Results, No Guesswork

Directigen™ Flu A+B

Directigen™ Flu A+B

The First Differentiated Test for Influenza A and Influenza B

Now there's a rapid, *in vitro* assay that detects *and differentiates* influenza A and B viral antigens—in less than fifteen minutes! BD Directigen™ Flu A+B is the *first* rapid assay on the market that distinguishes influenza A viral antigens from those of influenza B. Differentiation of influenza A and B can help physicians select appropriate, cost-effective treatment.

Directigen™ Flu A+B offers:

- * **A differentiated result**
- * **High accuracy**
- * **Fifteen minute time-to-results**
- * **Easy workflow**
- * **Ability to test a wide range of specimen types**



Directigen™ Flu A+B Rapid Test Kit

Various specimen types include nasopharyngeal wash, nasopharyngeal aspirate, nasopharyngeal swab, lower nasal swab, throat swab and bronchoalveolar lavage.

The Directigen™ Flu A+B rapid test is the latest addition to the BD line of rapid respiratory diagnostics that also includes RSV and Group A Streptococcus.



Directigen™ Flu A+B

20 Test Kit

Catalog No. 256010

Sensitivity and Specificity

Influenza A, Sensitivity of 95.7% and Specificity of 91.4% for nasopharyngeal aspirates.

Influenza B, Sensitivity of 87.5% and Specificity of 98.1% for nasopharyngeal aspirates.



BD Diagnostic Systems

7 Loveton Circle
Sparks, MD 21152-0999 USA
800.638.8663
www.bd.com/clinical

11 rue Aristide Bergès
38800 Le Pont de Claix, France
Tel: 33.4.7668.3636

2771 Bristol Circle
Oakville, Ontario
Canada L6H 6R5
Tel: 800.268.5430

30 Tuas Avenue 2
Singapore 639461
Tel: 65.6861.0633

Monte Pelvoux 111, 9th Floor
Col. Lomas de Chapultepec
11000 México D.F.
Tel: 52.5.237.1200

Rua Alexandre Dumas 1976
04717-004 São Paulo, S.P. Brazil
Tel: 55.11.5185.9833

Akasaka DS Building
5-26 Akasaka 8-chome
Minato-ku
Tokyo, 107 Japan
Tel: 81.24.593.5405