

Select Convenience  
and Quality for your  
Viral Collection and  
Transport Needs.

# BD *Cellmatics*<sup>TM</sup>

Viral Transport Pack



Indispensable to  
human health

# collection & transport

Maximize the quality of your virology lab results with the BD **Cellmatics** Viral Transport Pack. The BD **Cellmatics** Viral Transport Pack contains everything you need for viral collection and transport in a convenient package design.

Each BD **Cellmatics** Viral Transport Pack includes:

- > A screw-cap vial containing 2 mL **Difco**™ brand liquid viral transport medium that is tissue-culture compatible
- > Quality **Difco** brand viral transport medium composed of:
  - Hanks balanced salt solution with gelatin as a holding medium
  - HEPES biological buffer as used in cell culture media
  - Gentamicin sulfate to reduce the overgrowth of microorganisms
- > A cardboard transport pouch with printed instructions and patient ID label that protects the sample during transport



- > A sterile, rayon-tipped, plastic shaft swab that preserves the integrity of viruses without causing cell culture toxicity problems

**Cellmatics™** Viral Transport Pack  
50 packs . . . **252171**



## **BD Diagnostic Systems**

7 Loveton Circle  
Sparks, MD 21152-0999  
800.638.8663  
[www.bd.com/diagnostics](http://www.bd.com/diagnostics)

BD and BD Logo are trademarks of Becton, Dickinson and Company. Cellmatics and Difco are trademarks of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company. ©2001 BD 0-2528 October 2001 Printed in USA