



**BBL™ Sterile Pack
Plated Media**

**NOW FOR
ISOLATORS!**



Indispensable to
human health

BD BBL™ Isolator Pack Plated Media



Now the quality leader in plated media introduces new multi-wrap packaging that protects media from vaporized hydrogen peroxide exposure within an isolator.

Chemical Indicator Strip Assay for Vapor Hydrogen Exposure

RODAC™ Plates

POSITION OF CHEMICAL INDICATOR STRIP IN THE PACKAGING	OUTSIDE OF PACKAGING	INSIDE OF INNER BAG	INSIDE TOP PETRI DISH
NUMBER OF PACKAGE(S)	1	4	0
	2	4	0
	3	4	0
	4	4	0
	5	4	0
	6	4	0
	7	4	0
	8	4	0
	9	4	0
	10	4	0

100 mm Plates

POSITION OF CHEMICAL INDICATOR STRIP IN THE PACKAGING	OUTSIDE OF PACKAGING	INSIDE OF INNER BAG	INSIDE TOP PETRI DISH
NUMBER OF PACKAGE(S)	1	4	0
	2	4	0
	3	4	0
	4	4	0
	5	4	0
	6	4	0
	7	4	0
	8	4	0
	9	4	0
	10	4	0

4 represents complete VHP strip exposure with a color change (100% violet / 0% yellow)
0 represents no VHP strip exposure with no change in coloring (0% violet / 100% yellow)
Chemical indicator strip cat/lot/exp date: AMSCO Chemical Indicator CHEMDI-VHP™ Short Strip, Catalog No. NB302, Lot No. VHP-029901, Exp. Date 10/00

BD — Setting New Standards for the Industry

From the company that first developed irradiated prepared media — **BBL Isolator Pack Plated Media** was developed with the same quality standards you've come to depend on in the original **BBL Sterile Pack Plated Media**.

Validated Sterile

- Quarterly validations of **BBL Isolator Pack Plated Media** provide a high level of assurance against false positives
- Manufactured in an ISO 9000 Registered Quality System
- Produced in a cGMP manufacturing facility

Product/Description	Unit	Catalog No.
D/E Neutralizing Agar <i>Contains lecithin, polysorbate 80, sodium bisulfite to neutralize aldehydes, sodium thioglycollate to neutralize mercurials and sodium thiosulfate to neutralize chlorine.</i>	Package of 10.....	292645 Sterile RODAC Plates
	Bulk Pack of 100...	292646 Sterile RODAC Plates
	Package of 10.....	292647 Sterile Finger Dab™ Plates (100 mm x 15 mm style)
Trypticase™ Soy Agar with Lecithin and Polysorbate 80 <i>Contains lecithin to neutralize quaternary ammonium compounds and polysorbate 80 to neutralize phenolic disinfectants.</i>	Package of 10.....	299896 Sterile RODAC Plates
	Bulk Pack of 100...	292335 Sterile RODAC Plates
	Package of 10.....	292648 Sterile Finger Dab™ Plates (100 mm x 15 mm style)
	Bulk Pack of 100...	292649 Sterile Finger Dab™ Plates (100 mm x 15 mm style)
	Package of 5.....	292650 Sterile Finger Dab™ Plates (150 mm x 15 mm style)
	Trypticase™ Soy Agar <i>(Soybean-Casein Digest Agar)</i>	Package of 10.....
Bulk Pack of 100...		292652 Sterile Plates (100 mm x 15 mm style)
Package of 5.....		292272 Sterile Plates (150 mm x 15 mm style)
Sabouraud Dextrose Agar with Lecithin and Polysorbate 80	Package of 10.....	299937 Sterile RODAC Plates
	Bulk Pack of 100...	292337 Sterile RODAC Plates
	Package of 10.....	292653 Sterile Finger Dab™ Plates (100 mm x 15 mm style)
	Package of 5.....	292654 Sterile Finger Dab™ Plates (150 mm x 15 mm style)



Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152-9956 USA
800.638.8663
www.bd.com/biosciences/industrial

2771 Bristol Circle
Oakville, Ontario
Canada L6H 6R5
Tel: 800.268.5430

Monte Pelvoux 111, 9th Floor
Col. Lomas de Chapultepec
11000 México D.F.
Tel: 52.5.237.1200

Akasaka DS Building
5-26 Akasaka 8-chome
Minato-ku
Tokyo, 107 Japan
Tel: 81.24.593.5405

11 rue Aristide Bergès
38800 Le Pont de Claix, France
Tel: 33.4.7668.3636

30 Tuas Avenue 2
Singapore 639461
Tel: 65.8610.633

Rua Alexandre Dumas 1976
04717-004 São Paulo, S.P. Brazil
Tel: 55.11.5185.9833

All trademarks are the property of Becton, Dickinson and Company except CHEMDI-VHP, which is a trademark of Steris Corporation. ©2000. 1-2424 September 2000 Printed in USA