



Helping all people
live healthy lives

BD QTest™ Strep— Quick and Reliable Results in Your Office

*Only two minutes
from start to finish*



BD QTest™ Strep



Clearly, a better strep test.

The Fastest Test Available

Convenient Flexibility

Clear, Definitive Results


Simple Procedure

State-of-the-Art Accuracy

- ▲ Only two minutes from start to finish.
- ▲ Flow-through speed.
- ▲ Room temperature storage option.
- ▲ Samples may be stored and tested later.
- ▲ Procedure may be interrupted and continued at a later time.
- ▲ Distinct pink triangle clearly identifies positive test.
- ▲ (-) sign identifies negative test.



POSITIVE



NEGATIVE

- ▲ Extract sample (1 minute), add liposomes, add wash.
- ▲ Color-coded reagents.
- ▲ Built-in positive and negative procedural controls.
- ▲ Includes all necessary components for complete testing.
- ▲ Advanced liposome technology delivers accurate results.
- ▲ Positive and Negative Controls included in kit for convenience and quality assurance.

BD LINK2™

LINK2™ Strep A Test Kit
25 tests, Catalog No. 422012

BD QTest™

Available in 20, 40 and 80 test kits:

20 tests, Catalog No. 494775
40 tests, Catalog No. 494776
80 tests, Catalog No. 494780

For more information, contact your local BD distributor or call toll free 1-800-638-8663.



BD Diagnostics
7 Loveton Circle
Sparks, MD 21152-9956 USA
800.638.8663
www.bd.com/clinical

BD, BD Logo and all trademarks are the property of Becton, Dickinson and Company. ©2004 BD.

2-1003 January 2004 Printed in USA