

BD Sterile Pack Media Products

The largest breadth of line for your cleanroom needs



*Featuring the quality and dependability
of BBL™ and Difco™ media.*



BD Sterile Pack Media Products

Featuring unique packaging that minimizes the risk of false contamination both going into and coming out of your critical environment.

BD BBL™ Sterile Pack and Isolator Pack Prepared Plated Media

BD offers BBL™ Sterile Pack and BBL™ Isolator Pack Plated Media, specially designed for use in cleanrooms and isolators. A variety of specialized packaging is available, including standard single wrap; gamma-irradiated double wrap; and gamma-irradiated triple wrap which is resistant to vapor hydrogen peroxide. Complete validation packages available. The RODAC™ plate for surface sampling features a loose fitting lid that provides fast, easy one-hand manipulation. The RODAC™ SL is identical except with the addition of three lugs on the base providing a tight fit thereby reducing accidental contamination.



BD Sterile Pack Swabs

BD Sterile Pack Swabs are the first, double-wrapped, gamma irradiated, ready-to-use sterile swabs for surface sampling, combined with a rinse solution-filled tube.

BD Sterile Pack Bottles

BD Sterile Pack Bottles are terminally sterilized inside of autoclavable double-bags, resulting in a bottle exterior that is free from environmental contaminants and particulate matter.

	Units	Cat. Number	
		RODAC™	RODAC™ SL
BD BBL™ Sterile Pack Plated Media			
Trypticase™ Soy Agar with Lecithin and Polysorbate 80	10	221238	222248
Trypticase™ Soy Agar with Lecithin and Polysorbate 80	100	222207	222247
Trypticase™ Soy Agar with Lecithin, Polysorbate 80 and Penicillinase	10	221234	222246
D/E Neutralizing Agar	10	221232	222243
D/E Neutralizing Agar	100	222209	222251
D/E Neutralizing Agar with Penicillinase	10	292227	
Sabouraud Dextrose Agar	10	221235	
Sabouraud Dextrose Agar with Lecithin and Polysorbate 80	10	221233	222244
CDC Anaerobic Blood Agar	10	292229	
Rose Bengal Agar with Penicillinase	10	292225	

BD Sterile Pack Swabs

10 mL Sterile Rinse Solution	200	220518
------------------------------	-----	--------

BD Sterile Pack Bottles

100 mL Tryptic Soy Broth	10	257213
100 mL Fluid Thioglycollate Medium	10	257217



BD Diagnostics

7 Loveton Circle
Sparks, MD 21152-0999 USA
800.638.8663
www.bd.com/ds

2771 Bristol Circle
Oakville, Ontario
Canada L6H 6R5
Tel: 800.268.5430

Monte Pelvoux 111, 9th Floor
Col. Lomas de Chapultepec
11000 México D.F.
Tel: 52.55.59.99.82.00

Akasaka DS Building
5-26 Akasaka 8-chome
Minato-ku
Tokyo, 107 Japan
Tel: 81.24.593.5405

11 rue Aristide Bergès
38800 Le Pont de Claix, France
Tel: 33.4.7668.3636

30 Tuas Avenue 2
Singapore 639461
Tel: 65.6861.0633

Rua Alexandre Dumas 1976
04717-004 São Paulo, S.P. Brazil
Tel: 55.11.5185.9833