

# New BD™ Wide Mouth Sterility Bottles

## BD Media Solutions for Sterility Testing



W I D E  
M O U T H



Helping all people  
live healthy lives

# New BD™ Wide Mouth Sterility Bottles

- Sterility Assurance Level of 10<sup>-6</sup> reduces risk of false positives, saving time and money
- Terminally autoclaved
- No product contact with corrugate
- 48 mm opening permits better access for testing of articles or filters

- Conforms to Harmonized USP<sup>1</sup>, EP<sup>2</sup> and JP<sup>3</sup> requirements for sterility tests and to CLSI<sup>4</sup> requirements for recovery
- Clear FTM conforms to pharmacopeial formulation requirements and makes results easier to interpret

Cat #	Fill Volume	Format	Specifications		Qty / Case	SAL
			USP/EP/JP	CLSI		
<b>Tryptic Soy Broth</b>						
257486	100 mL	Wide Mouth Bottle	X	X	10	10 <sup>-6</sup>
257488	90 mL	Wide Mouth Bottle	X	X	10	10 <sup>-6</sup>
<b>Fluid Thioglycollate Medium</b>						
257485	100 mL	Wide Mouth Bottle	X	X	10	10 <sup>-6</sup>

Bottle Dimensions				
Capacity	Height	Outer Diameter	Inner Diameter	Cap
150 mL	130 mm	56 mm	48 mm	Metal twist-off

## Microbiology – It's what we do.

### Find out what we can do for you.

For more information on BD media solutions or to place an order, please contact your local BD sales manager or visit us on the web at [www.bd.com/ds](http://www.bd.com/ds).



Wide Mouth

1 United States Pharmacopeia Convention, Inc. 2009. The United States Pharmacopeia 32/The National Formulary 27, 2009. The United States Pharmacopeial Convention, Rockville, MD.  
2 European Pharmacopeia, 6th Edition. European Directorate for the Quality of Medicine, Council of Europe, 226 Avenue de Colmar BP907, F-67029 Strasbourg Cedex 1, France.  
3 Japanese Pharmacopeia, Fifteenth Edition.  
4 Clinical and Laboratory Standards Institute. 2004. Quality control for commercially prepared microbiological culture media; Approved Standard – Third Edition, M22-A3. Clinical and Laboratory Standards Institute, Wayne, PA.

