

BD ProbeTec™ ET Specimen Collection and Wet Transport System

Reduces Hands-on Time Up To 30%

The BD ProbeTec™ ET Collection Kit for CT/GC offers a unique design where specimens arrive ready for processing, eliminating swab transfer and manual expression of swab.

The BD Wet Swab Transport System combines: a pre-scored swab, pre-filled diluent tube and a specially designed cap which captures the swab and provides a leak-proof seal for convenient transportation to the test site while increasing user safety.

A System That Saves Time and Labor

- Specimen extraction performed in transport tube, eliminates additional steps
- Reduces number of capping/uncapping steps, saves time
- Increases testing flexibility with 4-6 day transport and storage at 2-27°C or up to 30 days storage at 2-8°C
- Color-coded system eases identification, reduces repeat specimen collection
- Large cleaning swab easily differentiated from collection swab, reduces specimen rejections

BD ProbeTec™ ET Wet Transport System

Pink Package – BD ProbeTec™ ET CT/GC Collection Kit for Endocervical Specimens **220142**

Blue Package – BD ProbeTec™ ET CT/GC Collection Kit for Urethral Specimens. . **220143**



Helping all people live healthy lives

BD ProbeTec™ ET Specimen Collection and Dry Transport System

A Convenient Choice for Your Clinician

For collection and transport of CT/GC, the BD ProbeTec™ ET Specimen Collection and Dry Transport System offers a simple-to-use system for the clinician. Collect the specimen, place swab in transport tube and send to lab, no additional steps necessary.

- Increases testing flexibility with 4-6 day transport and storage at 2-27°C or up to 30 days storage at 2-8°C
- Security friction cap enhances user safety by preventing sample “splashes”
- Color-coded system eases identification, reduces repeat specimen collection
- Large cleaning swab easily differentiated from collection swab, reduces specimen rejections

BD ProbeTec™ ET Dry Transport System

Pink Cap – BD ProbeTec™ ET CT/GC Endocervical Specimen Collection and Dry Transport System **440476**

Blue Cap – BD ProbeTec™ ET CT/GC Urethral Specimen Collection and Dry Transport System **440461**



BD Diagnostics
7 Loveton Circle
Sparks, MD 21152-0999
800.638.8663
www.bd.com/ds