

● Liquid-based cytology

● Computer-guided screening

● Immuno-chemistry



Helping all people
live healthy lives



Trust BD SurePath™

The right material for the right start

Only 3 easy steps to take a smear

1. Collect



Collect the cytology sample using either a broom-like device or a combination brush/plastic spatula with detachable heads.

2. Prepare



Drop the detachable head into the BD SurePath™ vial. Place the cap on the vial and tighten.


3. Send




Label and send the BD SurePath™ vial to the lab for processing.

The collection devices, recommended to be used with BD SurePath™, have all a detachable head, making it simple to collect, prepare and send a complete sample for analysis.

- You can be sure your samples **have been collected in a standardized way**
- You can be sure **that 100% of the collected cells are sent to the laboratory**. In contrast, using a rinse and swirl technique can result in an average 37% loss of cells ⁽¹⁾.



- You can be sure to **obtain a readable result** from your samples including bloody and inflammatory samples, as well as samples from menopausal patients or patients on birth control ⁽²⁾.



With BD SurePath™, you don't need to worry about discarding cells that could mean the difference between finding dysplasia and missing it! ⁽¹⁾.

References:

- (1) Bigras, G et al. J of Lower Genital Tract Disease. Vol 7 (3), 2003: 168-174.
 (2) Sweeney BJ, et al. Cancer Cytopathology. Vol 108 (1), 2006: 27-31.



BD Diagnostics TriPath

Erembodegem-dorp 86
 9320 Erembodegem, Belgium
 For Europe, Middle East and Africa
 Phone: +32 53 720 211
 Email: TriPathEurope@BD.com
 www.bd.com/TriPath