

Tubes BD Vacutainer^{MD} Plus

...la gamme de tubes en plastique BD Vacutainer^{MD}



Tubes citratés BD Vacutainer^{MD} Plus

BD sait que le remplissage des tubes citratés avec le volume d'échantillon recommandé est l'un des nombreux facteurs importants qui contribuent à l'obtention de résultats de qualité.

Vous avez demandé la présence d'un indicateur de remplissage pour vous aider à déterminer que le volume prélevé est suffisant.

Les tubes citratés BD Vacutainer^{MD} Plus portent maintenant un **indicateur de remplissage minimum** qui est moulé sur le tube. Cet indicateur correspond au volume minimum de sang requis pour effectuer une analyse de façon appropriée.

L'indicateur de remplissage nominal présent sur nos autres tubes BD Vacutainer^{MD} Plus demeure inchangé.

Les tubes sont destinés à être remplis jusqu'à l'indicateur de remplissage environ. Il est important de noter que la position de l'étiquette peut varier légèrement d'un tube à l'autre. Par conséquent, les indicateurs de remplissage ne montrent qu'approximativement le volume à prélever.



EDTA



EDTA



Glucose




Héparine



Plomb



PST^{MD}



Indicateur de remplissage nominal



Pour que chacun
vive en bonne santé

Courriel:
CustSvcCan@bd.com

Ou visitez-nous en ligne à
www.bd.com/ca

Services techniques BD
1.800.268.5430 poste 5567

BD Vacutainer® Plus Tubes

...the plastic product array of BD Vacutainer® Tubes



BD Vacutainer® Plus Citrate Tubes

BD knows that filling citrate tubes with the recommended sample volume is one of the many important factors in obtaining quality results.

You have requested the presence of a fill indicator to help you determine if the volume collected is sufficient.

The Citrate Plus Tubes now have a **minimum fill indicator** molded on the tube. This indicator represents the minimum volume of blood required for appropriate analysis.



Helping all people live healthy lives

The nominal fill indicator on our other BD Vacutainer® Plus Tubes remains unchanged.

Tubes are intended to fill approximately to the fill indicator. It is important to note that label placement can vary slightly from tube to tube. For this reason, fill indicators are approximate indicators of expected fill volume.



EDTA



EDTA



Glucose




Heparin



Lead



PST™



e-mail:

CustSvcCan@bd.com

Or visit us anytime online at
www.bd.com/ca

BD Technical Services at
1.800.268.5430 ext. 5567