

Product Description

Ref	Number	Needle Size	Description	Pkg. Qty.
-----	--------	-------------	-------------	-----------

BD Lo-Dose™ 3/10 mL Insulin Syringe with Permanently Attached Needle

SELF-CONTAINED				
328438	31 G x 5/16 in.	U-100 BD Ultra-Fine™ Short		100/Shelf Pack
328440	31 G x 5/16 in.	U-100 BD Ultra-Fine Short 1/2 Unit Scale		500/Case
309301	29 G x 1/2 in.	U-100 BD Ultra-Fine		200/Shelf Pack
309300	28 G x 1/2 in.	U-100 BD Micro-Fine™ IV (Orange)		1600/Case

BD Lo-Dose™ 1/2 mL Insulin Syringe with Permanently Attached Needle

SELF-CONTAINED				
328468	31 G x 5/16 in.	U-100 BD Ultra-Fine Short		100/Shelf Pack
329465	29 G x 1/2 in.	U-100 BD Ultra-Fine (Orange)		500/Case
309306	29 G x 1/2 in.	U-100 BD Ultra-Fine		200/Shelf Pack
309305	28 G x 1/2 in.	U-100 BD Micro-Fine IV (Orange)		1600/Case
BLISTER PACKAGE				
329461	28 G x 1/2 in.	U-100 BD Micro-Fine IV (Orange)		100/Shelf Pack 500/Case

BD™ 1 mL Insulin Syringe with Permanently Attached Needle

SELF-CONTAINED				
328418	31 G x 5/16 in.	U-100 BD Ultra-Fine Short		100/Shelf Pack
329410	28 G x 1/2 in.	U-100 BD Ultra-Fine (Orange)		500/Case
329412	27 G x 5/8 in.	U-100 BD Micro-Fine IV (Orange)		
309311	29 G x 1/2 in.	U-100 BD Ultra-Fine		200/Shelf Pack
309309	28 G x 1/2 in.	U-100 BD Micro-Fine IV (Orange)		1600/Case
309310	27 G x 5/8 in.	U-100 BD Micro-Fine IV (Orange)		
BLISTER PACKAGE				
329424	28 G x 1/2 in.	U-100 BD Micro-Fine IV (Orange)		100/Shelf Pack
329420	28 G x 1/2 in.	U-100 BD Micro-Fine IV Single Unit Graduation		500/Case

BD™ 1 mL Insulin Syringe with Detachable Needle

BLISTER PACKAGE				
329652	26 G x 1/2 in.	U-100 (Orange/Tan)		100/Shelf Pack
329622	25 G x 1 in.	U-100 (Orange/Blue)		800/Case
329651	25 G x 5/8 in.	U-100 (Orange/Blue)		

BD™ 1 mL Insulin Syringe Only

329650		U-100 Slip Tip		100/Shelf Pack
309629		U-100 BD Luer-Lok™ Tip		800/Case



BD Lo-Dose
3/10 mL Insulin Syringe



BD Lo-Dose
1/2 mL Insulin Syringe



BD 1 mL Insulin Syringe with
Permanently Attached Needle

1/2 mL Syringes
 1 mL Syringes
 3 mL Syringes

Section F
 Hypodermic Products/Diabetes

Product Description

Ref Number	Needle Size	Description	Pkg. Qty.
BD™ 1/2 mL Tuberculin Syringe with Permanently Attached Needle*			
305620	27 G x 1/2 in.	Regular Bevel	100/Shelf Pack 500/Case
BD™ 1 mL Tuberculin Syringe with Detachable Needle, Slip Tip*			
309623	27 G x 1/2 in.		100/Shelf Pack
309625	26 G x 3/8 in.	Intradermal Bevel	800/Case
309626	25 G x 5/8 in.		
309624	21 G x 1 in.		
BD™ 1 mL Tuberculin Syringe Only*			
309602	Slip Tip		100/Shelf Pack 800/Case
BD™ 1 mL Bulk Sterile Convenience Tray Pak*			
309701	Luer Slip Tip	(25 Syringes only per sterile tray, 12 trays per case of 300)	300/Case
BD™ 1 mL Allergy Syringe with Permanently Attached Needle*			
305500	28 G x 1/2 in.	Regular Bevel	100/Shelf Pack 500/Case
BD™ 1/2 mL Allergist Tray with Permanently Attached Needle* (25 Units Per Tray)			
305535	27 G x 1/2 in.	Regular Bevel	1000/Case
305536	27 G x 3/8 in.	Intradermal Bevel	
BD™ 1 mL Allergist Tray with Permanently Attached Needle* (25 Units Per Tray)			
305537	26 G x 1/2 in.	Regular Bevel	1000/Case
305538	26 G x 1/2 in.	Intradermal Bevel	
305539	26 G x 3/8 in.	Intradermal Bevel	
305540	27 G x 1/2 in.	Regular Bevel	
305541	27 G x 3/8 in.	Intradermal Bevel	
305542	27 G x 3/8 in.	Regular Bevel	
BD™ 1 mL Syringe/Needle Combination, Slip Tip*			
309597	26 G x 5/8 in.	1/100 mL Sub-Q	100/Shelf Pack 800/Case
BD Luer-Lok™ 1 mL Syringe Only			
309628	BD Luer-Lok 1/100 mL		100/Shelf Pack 800/Case
BD Luer-Lok™ 1 mL Syringe/Needle Combination			
309637	20 G x 1 in.		100/Shelf Pack 800/Case
BD™ 3 mL Syringe Only*			
309585	BD Luer-Lok Tip 1/10 mL		100/Shelf Pack
309586	Slip Tip		800/Case



BD 1 mL
 Tuberculin Syringe with Needle



BD 1 mL Allergist Tray



BD 1 mL
 Luer-Lok Syringe



BD 3 mL
 Luer-Lok Syringe

• New Product * Contains No Natural Rubber Latex

Product Description

Ref Number	Needle Size	Description	Pkg. Qty.
------------	-------------	-------------	-----------

BD™ 3 mL Syringe/Needle Combination BD Luer-Lok™ Tip*

309587	26 G x 5/8 in.	1/10 mL Sub-Q Thin Wall	100/Shelf Pack
309570	25 G x 5/8 in.	1/10 mL	800/Case
309581	25 G x 1 in.	1/10 mL	
309582	25 G x 1 1/2 in.	1/10 mL	
309588	23 G x 1 in.	1/10 mL IM Thin Wall	
309571	23 G x 1 in.	1/10 mL	
309589	23 G x 1 1/2 in.	1/10 mL IM Thin Wall	
309569	22 G x 3/4 in.	1/10 mL	
309572	22 G x 1 in.	1/10 mL	
309573	22 G x 1 1/4 in.	1/10 mL	
309574	22 G x 1 1/2 in.	1/10 mL	
309598	21 G x 1 in.	1/10 mL IV Thin Wall	
309575	21 G x 1 in.	1/10 mL	
309599	21 G x 1 1/2 in.	1/10 mL IV Thin Wall	
309577	21 G x 1 1/2 in.	1/10 mL	
309578	20 G x 1 in.	1/10 mL	
309579	20 G x 1 1/2 in.	1/10 mL	
309580	18 G x 1 1/2 in.	1/10 mL	

BD™ 3 mL Syringe With Blunt Fill Needle*

305060	BD Luer-Lok Tip 1/10 mL	100/Shelf Pack
	18 G x 1 1/2 in. Blunt Fill Needle	800/Case

BD™ 3 mL Bulk Sterile Convenience Tray Pak*

309702	BD Luer-Lok Tip (25 Syringes only per sterile tray, 12 trays per case of 300)	300/Case
--------	--	----------

BD™ 5 mL Syringe Only*

309603	BD Luer-Lok Tip 1/5 mL	100/Shelf Pack
301603	Slip Tip	400/Case

BD™ 5 mL Syringe/Needle Combination BD Luer-Lok™ Tip*

309630	22 G x 1 in.	1/5 mL	100/Shelf Pack
309631	22 G x 1 1/2 in.	1/5 mL	400/Case
309632	21 G x 1 in.	1/5 mL	
309633	21 G x 1 1/2 in.	1/5 mL	
309634	20 G x 1 in.	1/5 mL	
309635	20 G x 1 1/2 in.	1/5 mL	


BD™ 5 mL Syringe With Blunt Fill Needle*

305062	BD Luer-Lok Tip 1/5 mL	100/Shelf Pack
	18 G x 1 1/2 in. Blunt Fill Needle	400/Case

BD™ 5 mL Bulk Sterile Convenience Tray Pak*

309703	BD Luer-Lok Tip (25 Syringes only per sterile tray, 12 trays per case of 300)	300/Case
--------	--	----------

Types of Syringe Tips



BD 5 mL Syringe/Needle
Combination
BD Luer-Lok Tip

BD 5 mL Syringe

Product Description

Ref Number	Needle Size	Description	Pkg. Qty.
------------	-------------	-------------	-----------

BD™ 10 mL Syringe Only*

309604	BD Luer-Lok™ Tip 1/5 mL		100/Shelf Pack
301604	Slip Tip		400/Case
305482	Eccentric Tip		
309695	BD Luer-Lok Tip Control Syringe		25/Shelf Pack 100/Case
305462	Lab Syringe Eccentric Tip Non-Sterile, Clean, Ready-to-use, Bulk-Pak		100/Shelf Pack 400/Case

BD™ 10 mL Syringe/Needle Combination BD Luer-Lok™ Tip*

309640	22 G x 1 in.	1/5 mL	100/Shelf Pack
309642	21 G x 1 in.	1/5 mL	400/Case
309639	21 G x 1 1/2 in.	1/5 mL IV Thin Wall	
309643	21 G x 1 1/2 in.	1/5 mL	
309644	20 G x 1 in.	1/5 mL	
309645	20 G x 1 1/2 in.	1/5 mL	

BD™ 10 mL Syringe with Blunt Fill Needle*

305064	BD Luer-Lok Tip 1/5 mL 18 G x 1 1/2 in. Blunt Fill Needle		100/Shelf Pack 400/Case
--------	--	--	----------------------------

BD™ 10 mL Syringe Bulk Sterile Convenience Pak*

309605	BD Luer-Lok Tip 1 mL, 3/4 oz. in 1/8 oz. (20 Syringes only per sterile tray, 12 trays per case of 240)		240/Case
--------	---	--	----------

BD™ 20 mL Syringe Only*

309661	BD Luer-Lok Tip 1 mL, 3/4 oz. in 1/8 oz.		40/Shelf Pack
301625	Slip Tip 1 mL, 3/4 oz. in 1/8 oz.		160/Case
300613	Eccentric Tip		120/Shelf Pack 480/Case

BD™ 20 mL Syringe Bulk Sterile Convenience Pak Tray*

305617	BD Luer-Lok Tip (10 Syringes only per sterile tray, 12 trays per case of 120)		120/Case
--------	--	--	----------

BD™ 30mL Syringe Only*

309650	BD Luer-Lok Tip 1 mL, 1 oz. in 1/8 oz.		40/Shelf pack
309651	Slip Tip 1 mL, 1oz. in 1/8 oz.		160/Case

BD™ 30 mL Syringe Bulk Sterile Convenience Pak Tray*

305618	BD Luer-Lok Tip (10 Syringes only per sterile tray, 12 trays per case of 120)		120/Case
--------	--	--	----------



BD 10 mL Syringe/Needle Combination
BD 10 mL Syringe
BD Luer-Lok Tip



BD 10 mL Syringe
Bulk Sterile
Convenience Pak



BD 20 mL Syringe BD 30 mL Syringe

Product Description

Ref Number	Needle Size Description	Pkg. Qty.
------------	-------------------------	-----------

BD™ 60 mL Syringe Only*

309653	BD Luer-Lok™ Tip 1 mL, 2 oz. in 1/4 oz.	40/Shelf Pack
309654	Slip Tip 1 mL, 2 oz. in 1/4 oz.	160/Case
300866	Eccentric Tip 1 mL, 2 oz. in 1/4 oz.	60/Shelf Pack 240/Case

BD™ 60 mL Syringe Bulk Sterile Convenience Pak Tray*

309680	BD Luer-Lok Tip (20 Syringes only per sterile tray, 6 trays per case of 120)	120/Case
--------	---	----------

BD™ 2 oz. Catheter Tip Syringe Only*

309620	Catheter Tip 1 mL, 2 oz. in 1/4 oz.	40/Shelf Pack 160/Case
--------	-------------------------------------	---------------------------

BD™ Bulk Syringes, Non-Sterile, Clean, Ready-to-use*

Assembled Syringe with Scale, Siliconized, with BD Logo.
Bulk packed in double wrapped poly bags, in case carton.
Can be autoclaved or gas sterilized out of poly bags.
NOTE: Before autoclaving, remove plunger rod from barrel.

301025	1 mL Slip Tip	3000/Case
301073	3 mL BD Luer-Lok Tip	1600/Case
301077	3 mL Slip Tip	
301027	5 mL BD Luer-Lok Tip	1400/Case
301028	5 mL Slip Tip	
301029	10 mL BD Luer-Lok Tip	850/Case
301030	10 mL Slip Tip	
301031	20 mL BD Luer-Lok Tip	325/Case
301032	20 mL Slip Tip	
301033	30 mL BD Luer-Lok Tip	225/Case
301034	30 mL Slip Tip	
301035	60 mL BD Luer-Lok Tip	125/Case
301036	60 mL Slip Tip	
301037	2 oz. Catheter Tip	


BD™ Blunt Fill Needle*

305180	18 G x 1 1/2 in. Blunt Fill Needle	100/Shelf Pack 1000/Case
--------	------------------------------------	-----------------------------

BD™ Specialty Use Needles*

305106	30 G x 1/2 in.	100/Shelf Pack
305128	30 G x 1 in.	1000/Case
305136	27 G x 1 1/4 in.	
305124	25 G x 7/8 in.	
305120	23 G x 1 1/4 in.	
305197	16 G x 1 in.	
305198	16 G x 1 1/2 in.	
305129	21 G x 2 in.	

Types of Syringe Tips



Blunt Fill Needle

Product Description

Ref Number	Needle Size	Description	Pkg. Qty.
------------	-------------	-------------	-----------


BD PrecisionGlide™ Hypodermic Needle Tubing Specifications

Gauge	O.D.	Regular I.D.	Thin Wall I.D.
30	0.0120"	0.0060"	0.0070"
29	0.0130"	0.0070"	0.0080"
28	0.0140"	0.0070"	0.0080"
27	0.0160"	0.0075"	0.0100"
26	0.0180"	0.0095"	0.0120"
25	0.0200"	0.0095"	0.0120"
24	0.0220"	0.0115"	0.0140"
23	0.0250"	0.0125"	0.0150"
22	0.0280"	0.0155"	0.0180"
21	0.0320"	0.0195"	0.0220"
20	0.0350"	0.0230"	0.0255"
19	0.0420"	0.0270"	0.0315"
18	0.0500"	0.0330"	0.0380"
17	0.0590"	0.0410"	0.0460"
16	0.0650"	0.0470"	0.0520"

Hub is Polypropylene Needle Cannula is Stainless Steel

BD Brand Single-Use Regular Needles

Needle cannula is made of stainless steel. Polypropylene hub is designed to fit any slip tip, reusable or single-use. Strips of five sterile, fiber-free peel-apart packages. Shield protecting point and cannula. Package is color coded according to needle gauge. Lot control number on package.



BD™ Regular Bevel*

305109	27 G x 1/2 in.		100/Shelf Pack
305110	26 G x 3/8 in.	Intradermal Bevel	1000/Case
305111	26 G x 1/2 in.		
305122	25 G x 5/8 in.		
305125	25 G x 1 in.		
305127	25 G x 1 1/2 in.		
305143	23 G x 3/4 in.		
305145	23 G x 1 in.		
305155	22 G x 1 in.		
305156	22 G x 1 1/2 in.		
305165	21 G x 1 in.		
305167	21 G x 1 1/2 in.		
305115	26 G x 5/8 in.	Sub-Q	
305193	23 G x 1 in.	IM Thin Wall	
305194	23 G x 1 1/2 in.	IM Thin Wall	
305190	21 G x 1 1/2 in.	IV Thin Wall	
305177	21 G x 1 in.	IV Thin Wall	
305175	20 G x 1 in.		
305176	20 G x 1 1/2 in.		
305186	19 G x 1 in.	Thin Wall	
305187	19 G x 1 1/2 in.	Thin Wall	
305195	18 G x 1 in.		
305196	18 G x 1 1/2 in.		
305185	18 G x 1 1/2 in.	FILL Thin Wall	

Product Description

Ref Number	Needle Size	Description	Pkg. Qty.
------------	-------------	-------------	-----------

BD™ Short Bevel*

305159	22 G x 1 1/2 in.		100/Shelf Pack
305178	20 G x 1 in.		1000/Case
305179	20 G x 1 1/2 in.		
305188	19 G x 1 in.	Thin Wall	
305189	19 G x 1 1/2 in.	Thin Wall	
305199	18 G x 1 1/2 in.		


BD™ Bulk, Non-Sterile, Clean, Ready-to-use Needles with Shields*

Bulk packed 5000 units to a case in two (2) poly bags of 2500 units each			
REGULAR BEVEL			
303012	25 G x 1 in.		5000/Case
303010	25 G x 5/8 in.		
303011	23 G x 1 in.		
303009	22 G x 1 1/2 in.		
303008	21 G x 1 1/2 in.		
303007	20 G x 1 1/2 in.		
303006	19 G x 1 1/2 in.	Thin Wall	
303005	18 G x 1 1/2 in.		
SHORT BEVEL			
303015	22 G x 1 1/2 in.		
303013	18 G x 1 1/2 in.		


BD™ Syringe Accessories*

BD™ Syringe Storage Bins			
305677	Bin 10 3/4 in. x 8 1/4 in. x 7 in.		6/Case
305678	Bin 7 3/8 in. x 4 1/8 in. x 3 in.		24/Case
305679	Bin 5 3/8 in. x 4 1/8 in. x 3 in.		
305680	Divider for #305677		12/Case
305681	Divider for #305678		
305682	Divider for #305679		
BD™ Sterile Tip Cap Single			
305819	Sterile Syringe Tip Cap		200/Shelf Pack 2000/Case

Types of Needle Walls



The Thin Wall needle wall design is thinner than a Regular Wall needle. The inner diameter is bigger without any increase in the outside diameter. For example, a 21 gauge Thin Wall needle has the same outside diameter as a 21 gauge Regular Wall, but has the same inner diameter as a 20 gauge needle. Thin Wall allows the same flow volume as a standard walled needle one gauge larger. Thin Wall is available in 18, 19, 21, 22, and 23 gauge sizes.



Syringe Storage Bins


Product Description		
Ref Number	Needle Size Description	Pkg. Qty.

BD™ Diabetes Care Accessories

328233	BD Magni-Guide™ Insulin Syringe Magnifier	12/Case
325773	BD Ultra-Fine™ Sterile, Single Use, 30 G Lancets (100's)	100/Shelf Pack 2400/Case
325772	BD Ultra-Fine Sterile, Single Use, 30 G Lancets (200's)	200/Shelf Pack 2400/Case
322057	BD Ultra-Fine Sterile, Single Use, 33 G Lancets (100's)	100/Shelf Pack 2400/Case
323487	BD™ Home Sharps Container, 1.4 Quart	12/Case
326895	BD™ Alcohol Swab	100/Shelf Pack 1200/Case
328235	BD Safe-Clip™ Needle Clipping Device	12/Case
323488	BD™ Sharps Disposal by Mail, 1.4 Quart	

BD™ Pen Needles

328203	29 G x 12.7 mm	100/Shelf Pack
320109	31 G x 8 mm	1200/Case
320119	31 G x 5 mm	



BD Pen Needles