


# BD Beaver™ Otolaryngology Products



Helping all people  
live healthy lives

ENT/Head & Neck Procedure Blades

Myringotomy Blades and Knives

Tympanoplasty/Stapedectomy Blades

BD™ Custom ENT Procedure Kits/Packs

# BD Beaver™ Otolaryngology Products

## ENT/Head & Neck Procedure Blades

BD Beaver ENT/Head & Neck blades are designed to help enhance control and operative efficiency in a wide range of procedures, from facial nerve dissection to tonsillectomy. Because the requirements of each product are so unique and specialized, each BD BEAVER blade is shaped, angled and honed to help achieve the desired clinical result.

**376300 BD Beaver™ Mini-Blade™** 6 per box  
sharp on two sides, pointed tip



**376400 BD Beaver Mini-Blade** 24 per box  
rounded tip, sharp on one side



**376500 BD Beaver Mini-Blade** 24 per box  
smaller alternative to #11 general blade



**376700 BD Beaver Mini-Blade** 24 per box  
smaller alternative to #15 general blade



**376900 BD Beaver Mini-Blade** 6 per box  
straight, sharp all around, 180° cutting surface  
double bevel blade



**376910 BD Beaver Mini-Blade** 6 per box  
angled 60°, sharp all around, 180° cutting surface  
double bevel blade



**376800 Tonsillectomy Blade** 24 per box  
sickle sharp



**375940 Needle Blade** 3 per box  
1.5 mm cutting edge



**375910 Needle Blade** 3 per box  
3.0 mm cutting edge



**375920 Needle Blade** 3 per box  
5.0 mm cutting edge



**375930 Needle Blade** 3 per box  
7.0 mm cutting edge



**374769 Knife** 3 per box  
straight, sharp all around, 180° cutting surface  
double bevel blade



## BD Beaver™ Handles

**BD Beaver Handles will accommodate all Beaver blades.**

1 per box

**374310** Round, Knurled, Stainless 10.0 cm L x .50 cm D



**371350** Round, Knurled, Stainless 13.0 cm L x .50 cm D



# BD Beaver™ Otolaryngology Products

## Myringotomy Blades and Knives

BD Beaver Myringotomy blades and knives are intended to offer superior precision and control in delicate, typically pediatric myringotomy procedures. A narrower shaft has been engineered to improve visibility of the surgical field for greater intraoperative control. Tips designed for right or left cutting, or two-directional cutting, provide the versatility needed for fine tympanic incisions.

**377120 Myringotomy Blade** 6 per box  
spear tip, offset 45°, narrow shaft



**377122 Myringotomy Knife** 3 per box  
spear tip, offset 45°, narrow shaft



**377100 Myringotomy Blade** 6 per box  
lance tip, downward cut, offset 45°



**377110 Myringotomy Blade** 6 per box  
lance tip, upward cut, offset 45°



**377121 Myringotomy Blade** 6 per box  
lance tip, downward cut, narrow shaft, offset 45°



**377124 Myringotomy Knife** 3 per box  
lance tip, downward cut, narrow shaft, offset 45°



## Tympanoplasty/Stapedectomy Blades

A range of BD Beaver blade designs, shapes and angles is offered for the special requirements of both tympanoplasty and stapedectomy, from general canal procedures to elevating the delicate tympanic membrane. A narrower shaft has been engineered to enhance visibility of the operative site.

**377200 Tympanoplasty Blade** 3 per box  
60° angled, sharp all around  
2.5 mm W, bevel down



**374451 Tympanoplasty Blade** 3 per box  
50° angled, sharp all around  
2.0 mm W, bevel down



**377210 Tympanoplasty Blade** 3 per box  
straight, sharp all around  
2.5 mm W



**377300 Sickle Edge Blade** 3 per box  
straight 84 mm L  
5.0 mm cutting surface



## BD™ Custom ENT Procedure Kits/Packs



Featuring BD Beaver™ knives and blades

BD™ Custom ENT Procedure Kits and Packs can be customized to the surgeon's needs.



### **BD Medical**

Ophthalmic Systems

### **US Customer Service**

411 Waverley Oaks Road  
Waltham, MA 02452  
Tel: 781.906.7900  
Fax: 781.906.7992  
[www.bd.com/ophthalmology](http://www.bd.com/ophthalmology)

### **UK Customer Service**

The Danby Building  
Edmund Halley Road  
Oxford Science Park  
Oxford OX4 4DQ  
England  
Tel: 44.1.865.781666  
Fax: 44.1.865.781627

### **Contact in other European Countries**

11 rue Aristide Bergès - BP4  
38800 Le Pont de Claix  
France  
Tel: 33.4.76.68.32.27  
Fax: 33.4.76.68.32.87