




Indispensable to
human health

A background image showing three scientists in white protective suits, hoods, masks, and gloves. They are working in a laboratory setting. The scientist on the left is holding a petri dish with a red agar plate. The scientist in the middle is holding a small vial and a pipette. The scientist on the right is holding a pipette tip. The background is a gradient of blue and orange.

BD Sterile Pack Media Products

*The largest breadth of line for your cleanroom
needs, featuring the quality and dependability
of BBL™ and Difco™ media.*

BD Sterile Pack Media Products

Featuring unique packaging that minimizes the risk of false contamination both going into and coming out of your critical environment.

BD BBL™ Sterile Pack and Isolator Pack Prepared Plated Media

BD offers BBL™ Sterile Pack and BBL™ Isolator Pack Plated Media, specially designed for use in cleanrooms and isolators. A variety of specialized packaging is available, including gamma-irradiated double and triple wrap; and gamma irradiated wrap which is resistant to vapor hydrogen peroxide. Complete validation packages available.

BD Sterile Pack Swabs

BD Sterile Pack Swabs are the first, double-wrapped, gamma irradiated, ready-to-use sterile swabs for surface sampling, combined with a rinse solution-filled tube.

BD Sterile Pack Bottles

BD Sterile Pack Bottles are terminally sterilized inside of autoclavable double-bags, resulting in a bottle exterior that is free from environmental contaminants and particulate matter.



Catalog Number Units

BD BBL™ Sterile Pack Plated Media RODACT™ Plates

Trypticase™ Soy Agar with Lecithin and Polysorbate 80	221238/222207	10/100
Trypticase™ Soy Agar with Lecithin, Polysorbate 80 and Penicillinase	221234	10
D/E Neutralizing Agar	221232/222209	10/100
D/E Neutralizing Agar with Penicillinase	292227	10
Sabouraud Dextrose Agar	221235	10
Sabouraud Dextrose Agar with Lecithin and Polysorbate 80	221233	10
CDC Anaerobic Blood Agar	292229	10
Rose Bengal Agar with Penicillinase	292225	10

BD Sterile Pack Swabs

10 mL Sterile Rinse Solution	220518	200
------------------------------	--------	-----

BD Sterile Pack Bottles

100 mL Tryptic Soy Broth	257213	10
100 mL Fluid Thioglycollate Medium	257217	10
300 mL Fluid A	257263	10
300 mL Fluid D	254978	10

Please refer to our product catalogue for a complete product portfolio.
Customized products available on request.



BD Diagnostic Systems
www.bd.com/industrial

Local Offices

Austria Am Concorde Park E1/7 2320 Schwechat Tel: 43.1.70.63.66.00	Eastern Europe/ Middle East, Africa Tullastrasse 8-12 69126 Heidelberg Tel: 49.6221.305.184	Germany Tullastrasse 8-12 69126 Heidelberg Tel: 49.6221.305.0	Italy Via Caldera 21 20153 Milano Tel: 39.02.48.24.01	Spain Camino de Valdeoliva s/n San Agustin de Guadalix 28750 Madrid Tel: 34.91.848.81.00	Switzerland Binningerstrasse 94 4123 Allschwil Tel: 41.61.485.22.22
Belgium POB 13, Denderstraat 24 9320 Erembodegem Tel: 32.53.720.211	France 11 rue Aristide Bergès 38800 Le Pont de Claix Tel: 33.4.76.68.36.36	Greece 5 Amfitheas Avenue 17122 Nea Smyrni Athens Tel: 30.210.940.77.41	Netherlands Postbus 757 2400 AT Alphen aan den Rijn Tel: 31.20.582.94.16	Sweden Instrumentvägen 19 PO Box 32054 12611 Stockholm Tel: 46.8.775.51.00	U.K. 21 Between Towns Road Cowley Oxford OX4 3LY Tel: 44.1865.781.677