

Products designed to keep you safe

Protecting yourself while caring for others

BD Microtainer® Contact-Activated Lancet



Catalog No.	NDC/HRI No.	Lancet Size
366592	08290-3665-92	30 G 1.5 mm



BD Ultra-Fine™ Pen Needles and Insulin Syringes are covered at the preferred co-pay by virtually all managed healthcare plans, including Medicare Part D.†

To receive mail about new diabetes products and special offers and promotions, call 1.888.BDCARES†, Monday through Friday, 8 am to 8 pm Eastern Time.

†NOT FOR MEDICAL OR EMERGENCY INFORMATION.

†Co-pays and preferred status vary by plan. Check with your plan for the savings associated with a preferred brand.

BD SafetyGlide™ Insulin Syringe with permanently attached needle



Catalog No.	No. NDC/HRI No.	Syringe/Needle Size
305930	08290-5930-01	1 mL 29 G x 1/2"
305932	08290-5932-01	1/2 mL 29 G x 1/2"
305934	08290-3059-34	1/2 mL 30 G x 5/16"
305935	08290-3059-35	3/10 mL 29 G x 1/2"
305937	08290-3059-37	3/10 mL 31 G x 5/16"

BD Safety-Lok™ Insulin Syringe with permanently attached needle



Catalog	No. NDC/HRI No.	Syringe/Needle Size
329464	08290-3294-64	1 ml 29 G x 1/2"

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2010 BD 10D01012315CB3R

BD
1 Becton Drive
Franklin Lakes, NJ 07417
www.BD.com/us/diabetes

Comfort is priceless.

BD diabetes products for all your needs.



Helping all people live healthy lives



For more information go to www.BD.com/us/diabetes

BD provides the products that help you manage diabetes...comfortably

For over a century, BD has been making products that help people live healthier, happier and more comfortable lives. During that time, BD has expanded and refined its product line to help people better manage their diabetes.

Length: Insulin needles come in different lengths. BD Original Needles are 12.7 mm (1/2") long, and BD Short Needles are 8 mm (5/16") long. The BD Nano Needle 4mm (3/16") and BD Mini Needle 5mm (5/32") are available only in pen needles.

Gauge: Insulin needles are available in many gauges (G), or thicknesses. The higher the gauge, the thinner the needle. For example, a 32 G needle is actually thinner and more comfortable than a 29 G needle.

BD Insulin Syringes	
BD Ultra-Fine™ Short	8 mm (5/16") x 31 G
BD Ultra-Fine™	12.7 mm (1/2") x 30 G
BD Micro-Fine™	12.7 mm (1/2") x 28 G

BD Pen Needles	
BD Ultra-Fine™ Nano	4 mm (5/32") x 32 G
BD Ultra-Fine™ Mini	5 mm (3/16") x 31 G
BD Ultra-Fine™ Short	8 mm (5/16") x 31 G
BD Ultra-Fine™ Original	12.7 mm (1/2") x 29 G

Other helpful BD Items

BD™ Alcohol Swabs

For pre-injection and pre-sampling cleansing



Catalog No.	NDC/HRI No.
326895	08290-3268-95

BD Magni-Guide™

Insulin syringe magnifier and needle guide to help patients with vision or dexterity problems



Catalog No.	NDC/HRI No.
328233	08290-3282-33

BD Lancets

BD Ultra-Fine™ 33 Gauge Lancets

- No lancet is thinner
- Fits most lancing devices



Quantity	Catalog No.	NDC/HRI No.
100	322057	08290-3220-57


BD Ultra-Fine™ 30 Gauge Lancets

- Ultra-thin for greater comfort
- Fits most lancing devices



Quantity	Catalog No.	NDC/HRI No.
100	325773	08290-3257-73
200	325772	08290-3257-72

Why is Proper Disposal Important?



Proper syringe disposal will help to:

- Store and safely dispose of used syringes and lancets,
- Protect trash collectors from accidental needlesticks,
- Prevent your used syringes from falling into the wrong hands,
- Protect the environment.

Some things are Single Use Only

Why Insulin Syringes and Pen Needles shouldn't be reused:



Sterility

If you reuse a needle, you should be aware that it is no longer sterile

Injection pain

Insulin needles are finely ground and lubricated. Reuse leads to tip damage and loss of lubricant which increases pain.

Dosage accuracy

Keeping a pen needle on your insulin pen may affect the accuracy of your dose. Talk to your healthcare professional for more details.



BD Disposal

BD™ Home Sharps Container

For safe disposal of used insulin syringes, pen needles and lancets



Catalog No.

NDC/HRI No.

323487

08290-3234-87

BD™ Sharps Disposal by Mail

Sharps-filled containers are mailed back for destruction. Postage pre-paid.



Catalog No.

NDC/HRI No.

323488

08290-3234-88

BD Safe-Clip™

Needle clipping and storage device safely holds up to 1,500 needles



Catalog No.

NDC/HRI No.

328235

08290-3282-35

With more than 80 years of experience in making the world's leading brand of insulin delivery products, BD insulin syringes and pen needles meet your needs for quality and comfort. The needles are designed for the ultimate in comfort and come in a variety of lengths and gauges.

Syringe Volume:

In addition to gauge and length, there is one more factor to consider when choosing a syringe: volume. Choose the smallest syringe that can hold your prescribed amount of insulin. It's easier to read the markings on the smaller syringes.


The three syringe sizes are:

Capacity	For doses up to:
3/10 mL (Half-Unit Scale)	30 units
3/10 mL	30 units
1/2 mL	50 units
1 mL	100 units

Why shorter is more comfortable

Most people prefer short needles. They are more comfortable and less threatening, so it's easier to start injecting insulin. Short needles also help prevent injections from going into muscle, which can result in pain and poor blood sugar control.

Be sure to check with your healthcare professional before using a shorter needle.



BD Nano 4mm and BD Mini 5mm only available in pen needles

BD Ultra-Fine™ Pen Needles
fit all diabetes pens and dosers
sold in the U.S.*

*As of April 2010

BD Insulin Syringes

BD Ultra-Fine™ Short Needle Insulin Syringe



8 mm (5/16") x 31 G

Capacity	Catalog No.	NDC/HRI No.
1 mL (100 unit)	328418	08290-3284-18
1/2 mL (50 unit)	328468	08290-3284-68
3/10 mL (30 unit)	328438	08290-3284-38

BD Ultra-Fine™ Needle Insulin Syringe



12.7 mm (1/2") x 30 G

Capacity	Catalog No.	NDC/HRI No.
1 mL (100 unit)	328411	08290-3284-11
1/2 mL (50 unit)	328466	08290-3284-66
3/10 mL (30 unit)	328431	08290-3284-31

BD Ultra-Fine™ Short Needle Insulin Syringe with half-unit-scale



8 mm (5/16") x 31 G

Capacity	Catalog No.	NDC/HRI No.
3/10 mL (30 unit)	328440	08290-3284-40

BD Micro-Fine™ Needle Insulin Syringe



12.7 mm (1/2")x 28 G

Capacity	Catalog No.	NDC/HRI No.
1 mL (100 unit)	328410	08290-3284-10
1/2 mL (50 unit)	328465	08290-3284-65
3/10 mL (30 unit)	328430	08290-3284-30

BD Pen Needles

BD Ultra-Fine™ Nano Pen Needle



4 mm (5/32") x 32 G

Catalog No.	NDC/HRI No.
320122	08290-3201-22

BD Ultra-Fine™ Mini Pen Needle



5 mm (3/16") x 31 G

Catalog No.	NDC/HRI No.
320119	08290-3201-19

BD Ultra-Fine™ Short Pen Needle



8 mm (5/16") x 31 G

Catalog No.	NDC/HRI No.
320109	08290-3201-09

BD Ultra-Fine™ Original Pen Needle



12.7 mm (1/2") x 29 G

Catalog No.	NDC/HRI No.
328203	08290-3282-03

Never carry a diabetes pen with the pen needle attached

A pen needle has openings on both ends that allow for a possible exchange of insulin and air. When you leave a needle attached to your pen, changes in the environment could impact your dosage.