

BD Nano just got better! Now with PentaPoint™ Comfort



Helping all people
live healthy lives

A patented advance in injection comfort

BD's patented 5-bevel needle tip creates a flatter, thinner surface to help penetrate the skin with significantly greater ease for a smoother and gentler injection.



Even
Better!



**BD Nano just got better!
So why recommend anything else?**

Upgrade your patients today!

www.BD.com/nano

The BD Ultra-Fine™ Nano Pen Needle

4 mm (5/32") x 32G

Catalog #: 320122

NDC/HRI #: 08290-3201-22



The BD Ultra-Fine™ Mini Pen Needle

5 mm (3/16") x 31G

Catalog #: 320119

NDC/HRI #: 08290-3201-19



The BD Ultra-Fine™ Short Pen Needle

8 mm (5/16") x 31G

Catalog #: 320109

NDC/HRI #: 08290-3201-09



The BD Ultra-Fine™ Original Pen Needle

12.7 mm (1/2") x 29G

Catalog #: 328203

NDC/HRI #: 08290-3282-03



Needles are actual size

Compatibility Chart

BD Ultra-Fine™ Pen Needles fit all pens sold in the U.S.*

Drug Manufacturer	Pen Name	Drug	Drug NDC Number	Place Where Drug Manufacturers Reference BD Needles	Reference to BD Ultra-Fine™ Pen Needles
Eli Lilly and Company	HumaPen® LUXURA® HD 	Humalog® 3mL	00002-7516-59	User Manual, Website	"HumaPen® LUXURA™ HD is suitable for use with Becton, Dickinson and Company pen needles."
	Humalog® KwikPen™ 	Humalog® Humalog® Mix 75/25® Humalog® Mix 50/50	00002-8799-59 00002-8797-59 00002-8798-59	User Manual, Brochures, Prescribing Card	"When prescribing KwikPen™ prefilled with the Humalog brand of insulins you will need to write a separate prescription for BD needles."
	Humulin® Pen 	Humulin® N Humulin® 70/30	00002-8730-59 00002-8770-59	User Manual, Website	"This device is suitable for use with Becton, Dickinson and Company's insulin pen needles."
Sanofi-Aventis	Apidra® SoloStar® 	Apidra®	0088-2502-05	Package Insert Patient Instructions for use Patient Toll Free Phone Number	"SoloStar® is compatible with all pen needles from Becton Dickinson and Company"
	Lantus® SoloStar® 	Lantus®	0088-2219-05	Package Insert User Manual, Device Page, Website	"BD Ultra-Fine™ needles to be used in conjunction with SoloStar® and OptiClik® are sold separately and are manufactured by BD." "BD Ultra-Fine™ needles are compatible with SoloStar®. These are sold separately and are manufactured by BD."
Amylin Pharmaceuticals	Byetta® 	Byetta® 10MCG Byetta® 5MCG	66780-0212-01 66780-0210-07		
Amylin Pharmaceuticals	SymlinPen® 	Symlin® 60 Symlin® 120	66780-0115-02 66780-0121-02		
Novo Nordisk, Inc.	Victoza® 	Victoza®	N/A		
	NovoLog® Mix FlexPen® 	NovoLog® Mix 70/30	N/A		
	NovoPen® Junior 	Novo Nordisk PenFill 3ml	N/A		
	Novo Pen® 3 	Novo Nordisk PenFill 3ml	N/A		
	Levemir® FlexPen® 	Levemir®	N/A		
Owen Mumford	AutoPen® 	N/A	N/A		

Are the pen needles you carry ISO-certified?

ISO stands for the International Organization for Standardization (ISO), an international organization that develops precise processes to create an industry-wide consensus. International Standardization is a multi-stakeholder process involving industry, science, academia, consumer and governments.

BD pen needles are tested and certified for connectivity and dose accuracy with Insulin Pens in accordance with ISO standards.

*As of February, 2012.

†Tested with adults of BMI 20-49.

The BD Nano is guaranteed to be compatible with all diabetes pens sold in the US.*

BD Ultra-Fine™ Pen Needles are recommended by Eli Lilly and Sanofi-Aventis.

BD, BD Logo, BD Ultra-Fine and PentaPoint are trademarks are property of Becton, Dickinson and Company. ©2012 BD.

All other brands are trademarks of their respective owners.

120881275PNC