

BD Microtainer® MAP Microtube for Automated Process

When Every Drop Counts



Helping all people
live healthy lives

When caring for your most fragile patients, keep your focus on the patient and let **BD Microtainer® MAP Microtube for Automated Process** take care of the sample.


At BD we understand nursing's vital role in the patient's experience. We are proud to equip you with innovative solutions intended to help you deliver the quality of care your fragile patients deserve.

BD Microtainer® MAP Microtube for Automated Process is designed to help you keep the focus on your patient by enabling:

- Patient identification with standard label to avoid potential labeling errors
- Reduction of repeat specimen collection caused by non-standardized process, transfer and labeling errors
- Larger tube size: easier to handle one piece design eliminates the need for assembly or discard of extra contaminated parts
- Your fragile vein and anemic patients can now provide smaller samples which can be processed automatically due to the tube's innovative design

BD Microtainer® MAP Microtube for Automated Process

Twist-assist closure, for easier cap removal, **snaps closed** and ensures **safer leak-free** transportation to the lab



Pierceable cap, compatible with instrument probes

One piece 13x75 mm tube with **integrated collector**


Three clearly **visible fill lines** to ensure proper collection volume

Number of inversions illustrated on the tube label

Full size standard label with **V-Notch™** feature provides a visual guide for proper placement of labeling



Packaging is recyclable



Ref #	Description	Box/Case Quantities
363706	BD Microtainer® MAP Microtube for Automated Processing: K2EDTA Additive 1.0 mg, 13x75 mm, Volume 250-500 µL	50/200

We value your time. This product's innovative design can help you to eliminate costly work-arounds such as extenders, syringes, transfers, blank tubes for labels, and help you to avoid potential costly mistakes.

BD Microtainer® MAP Microtube for Automated Process is one of the many tools designed by BD to increase the trust between you and your patients in order to maintain uncompromising standards of quality, integrity and safety.

To learn more about **BD products** and services contact us at:

BD Global Technical Services at **1.800.631.0174**
 BD Customer Service at **1.888.237.2762**
 or visit **www.bd.com/vacutainer/MAP**



BD Diagnostics
 Preanalytical Systems
 1 Becton Drive
 Franklin Lakes, NJ 07417
www.bd.com/vacutainer