

A Part of the
BD Vacutainer® System of
Specimen Collection Products



Ordering Information		
Ref.#	Description	Packaging
BD Vacutainer® Eclipse™ Blood Collection Needle with Pre-Attached Holder		
368650	With Pre-Attached Holder, 21 G x 1 1/4 in.	100/Case
368651	With Pre-Attached Holder, 22 G x 1 1/4 in.	100/Case
BD Vacutainer® Eclipse™ Blood Collection Needle		
368607	21 G x 1 1/4 in.	48/Box 480/Case
368608	22 G x 1 1/4 in.	48/Box 480/Case
BD Vacutainer® One Use Holder		
364815	One Use Holder	250/Bag 1000/Case
BD Vacutainer® Safety-Lok™ Blood Collection Sets		
367281	With Luer, 21 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367283	With Luer, 23 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367285	With Luer, 25 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367287	With Luer, 21 G 3/4 in., 7 in. Tubing	50/Box 200/Case
367292	With Luer, 23 G 3/4 in., 7 in. Tubing	50/Box 200/Case
367294	With Luer, 25 G 3/4 in., 7 in. Tubing	50/Box 200/Case
367296	Without Luer, 21 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367297	Without Luer, 23 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367298	Without Luer, 25 G 3/4 in., 12 in. Tubing	50/Box 200/Case
BD Vacutainer® Safety-Lok™ Blood Collection Sets with Pre-Attached Holder		
368652	With Pre-Attached Holder, 21 G 3/4 in., 12 in. Tubing	25/Box 200/Case
368653	With Pre-Attached Holder, 23 G 3/4 in., 12 in. Tubing	25/Box 200/Case
BD Vacutainer® Alcohol Swabs Isopropyl Alcohol, 70% v/v		
365016	150 swabs per dispenser	10/Box 40/Case

Ordering Information		
Ref.#	Description	Packaging
BD Vacutainer® Push Button Blood Collection Sets		
367338	With Luer, 21 G 3/4 in., 7 in. Tubing	50/Box 200/Case
367344	With Luer, 21 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367326	Without Luer, 21 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367336	With Luer, 23 G 3/4 in., 7 in. Tubing	50/Box 200/Case
367342	With Luer, 23 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367324	Without Luer, 23 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367335	With Luer, 25 G 3/4 in., 7 in. Tubing	50/Box 200/Case
367341	With Luer, 25 G 3/4 in., 12 in. Tubing	50/Box 200/Case
367323	Without Luer, 25 G 3/4 in., 12 in. Tubing	50/Box 200/Case
BD Vacutainer® Push Button Blood Collection Sets with Pre-Attached Holder		
367352	With Holder, 21 G 3/4 in., 12 in. Tubing	20/Box 100/Case
368656	With Holder, 23 G 3/4 in., 12 in. Tubing	20/Box 100/Case
BD Vacutainer® Luer Adapter		
367290	Multiple Sample Luer Adapter	100/Box 1000/Case
BD Vacutainer® Stretch Latex-Free Tourniquet		
367203	Stretch Latex-Free Tourniquet	25/Box 500/Case

For more information on BD Safety-Engineered products:

BD Global Technical Services: **1.800.631.0174**

BD Customer Service: **1.888.237.2762**

www.bd.com/vacutainer



BD Diagnostics
Preanalytical Systems
1 Becton Drive
Franklin Lakes, NJ 07417

BD Vacutainer® Eclipse™

Blood Collection Needle **with Pre-Attached Holder**



Helping all people
live healthy lives

convenience
compliance
commitment

Convenience

Because practicing safety should be simple

The BD Vacutainer® Eclipse™ Blood Collection Needle with Pre-Attached Holder is ready to use right out of the package, with no assembly required.

BD Vacutainer® Eclipse™ Blood Collection Needle with Pre-Attached Holder

Offers a simple, effective way to collect blood, while reducing accidental needlestick injuries

- INTUITIVE and EASY TO USE, so there is minimal change to your technique
- NEEDLE “Bevel Up” orientation toward the pink safety shield PROVIDES EASY IDENTIFICATION of bevel position
- SAFETY SHIELD is an “INTEGRAL” part of the needle; eliminates the need for secondary safety accessory device
- Immediate ONE-HANDED safety feature ACTIVATION reduces exposure time and doesn’t require a hard surface for activation
- SAFETY SHIELD design allows HANDS to REMAIN BEHIND the NEEDLE at all times
- Uses BD PrecisionGlide™ needle technology to ensure HIGHEST QUALITY and NEEDLE SHARPNESS
- Each holder is engraved with the words “DO NOT REUSE” and “SINGLE USE ONLY” to serve as a reminder for single-use holder compliance



BD Vacutainer® Safety-Lok™
Blood Collection Set with
Pre-Attached Holder
#368653

BD Vacutainer® Eclipse™
Blood Collection Needle
with Pre-Attached Holder
#368650

BD Vacutainer® Push Button
Blood Collection Set with
Pre-Attached Holder
#368656

Compliance

Because your focus should be on your patient

The BD Vacutainer® Eclipse™ Blood Collection Needle with Pre-Attached Holder:

- Is designed to help your institution comply with OSHA SHIB 10-15-03*, by reducing the risk of holder reuse.
- Provides assurance that the number of holders will always equal the number of Eclipse™ Blood Collection Needles.
- Each pre-attached holder is engraved with the words “DO NOT REUSE” and “SINGLE USE ONLY” to serve as a reminder for single-use holder compliance.



- 1 Pink safety shield
- 2 Needle is always bevel up toward the pink safety shield
Same BD PrecisionGlide™
needle quality and sharpness

- 3 Thumb pad for easy, one-handed activation
- 4 Safety shield is an “integral” part of the needle

- 5 Pre-load line
- 6 Engraved text to help drive compliance

Commitment

Because your product should be reliable

- The BD Vacutainer® Eclipse™ Blood Collection Needle with Pre-Attached Holder protects the healthcare worker from nonpatient-end needle exposure.
- The Eclipse™ safety shield reduces needlestick risk throughout the hospital, for both hospital personnel and patients.
- After use, the needle can be easily and completely covered by the protective shield that is securely locked in place.
- Improved safety protection from time of activation to final disposal may help reduce long-term costs of needlestick injuries.

*Disposal of Contaminated Needles and Blood Tube Holders Used for Phlebotomy