

BD Vacutainer® Plus UA Preservative Tube

The quick, safe way to minimize preanalytical errors in urine testing



Helping all people
live healthy lives

BD Vacutainer® Plus UA Preservative Tube

The **BD Vacutainer® Plus UA Preservative Tube** is designed to give you more time to test and process your urinalysis samples. With its unique formulation samples can be preserved for up to **72 hours without refrigeration**.

This allows you to maintain sample integrity and helps to minimize preanalytical errors and unnecessary cultures due to bacterial contamination or overgrowth.

Clinically proven with the Siemens Clinitek® Atlas®, Iris® iQ200 and Sysmex® UF-50 instruments.

Additional benefits include:

- Improved accuracy and workflow
- Easy mixing
- Enhanced instrument compatibility

The BD Vacutainer® Plus UA Preservative Tube is also mercury free*.



*In compliance with the memorandum of understanding between the EPA and American Hospital Association, calling for the elimination of all mercury waste in hospitals by the year 2005.

BD Vacutainer® Plus UA Preservative Tube


Now featuring a spray-coated additive for **easier mixing**

Providing you with **enhanced** instrument compatibility

Before



After



Results of **Improved!** BD Vacutainer® UA Preservative Tube¹

Ref #	Description	Box/Case Quantities
364957	Complete Kit: Sterile Screw-Cap Collection Cup with Integrated Transfer Device and 8.0 mL, 16 x 100 mm Plus Plastic Conical Tube with Preservative for Urinalysis and 4.0 mL, 13 x 75 mm Plus Plastic C&S Preservative Tube and Castile Soap Towelettes	50 Kits/Case
364946	Urinalysis Cup Kit: Sterile Screw-Cap Collection Cup with Integrated Transfer Device and 8.0 mL, 16 x 100 mm Plus Plastic Conical Tube with Preservative for Urinalysis	50 Kits/Case
364992	Bulk Tube: 8.0 mL, 16 x 100 mm Plus Plastic Conical Bottom Tube with Preservative for Urinalysis	100/1000
364943	Urinalysis Transfer Straw Kit: Transfer Straw and 8.0 mL, 16 x 100 mm Plus Plastic Conical Tube with Preservative for Urinalysis	50/200

For additional information on the **BD Vacutainer® Urine Collection Products** or any of the **BD Vacutainer®** Family of products and services please feel free to contact us via:

BD Global Technical Services at **1.800.631.0174**
or e-mail **vacutainer_techservices@bd.com**

BD Customer Service at **1.888.237.2762**
or visit us anytime online at **www.bd.com/vacutainer**



BD Diagnostics
Preanalytical Systems
1 Becton Drive
Franklin Lakes, NJ 07417
www.bd.com/vacutainer

¹ A Comparative Evaluation of Spray Coated Urinalysis Preservative Tubes with Current BD Vacutainer® Plus Urinalysis Preservative Tubes for Routine Urinalysis Testing on the Sysmex® UF-50 and Siemens Clinitek® Atlas Systems.