

BD[®] NRFit[™] offering

EMEA detail sheet

BD[®] spinal NRFit[™] introducers and needles

BD[®] spinal NRFit[™] needles, BD[®] spinal introducer NRFit[™] needles and BD[®] spinal NRFit[™] needle sets maintain the same design, materials of construction (except the coloring) and manufacturing processes as BD[®] Luer spinal Products, except for design changes to the needle hub to comply with ISO 80369-6. The analogous Luer and Bulk Non-Sterile (BNS) product codes are provided below for reference.

BD[®] spinal introducer NRFit[™] needles

Code NRFit SSU	Description	Metric dimensions (mm)	Wall	BD manufacturing site	Analogous BNS code	Analogous Luer code
BD [®] Spinal Introducer NRFit [™] Needles						
400919	22 G x 1 ¼ in. BD [®] Spinal Introducer NRFit [™] Needle	0.7 x 31.8	Regular	San Agustín	400916	-
400177	20 G x 1 ¼ in. BD [®] Spinal Introducer NRFit [™] Needle	0.9 x 31.8	Regular	San Agustín	400076	-

BD[®] Whitacre spinal NRFit[™] needles

Code NRFit SSU	Description	Metric dimensions (mm)	Wall	BD manufacturing site	Analogous BNS code	Analogous Luer code
BD [®] Whitacre Spinal NRFit [™] Needles						
400982	27 G x 3 ½ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.4 x 88.9	Regular	San Agustín	400968	402051
400222	27 G x 4 ¼ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.4 x 103.2	Regular	Juncos	400232	405111
400223	27 G x 4 ⅝ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.4 x 119.1	Regular	Juncos	400969	409443
400980	25 G x 3 ½ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.5 x 88.9	Thin	San Agustín	400966	402050
400220	25 G x 4 ¼ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.5 x 103.2	Thin	Juncos	400231	405110
400221	25 G x 4 ⅝ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.5 x 119.1	Regular	Juncos	400967	409442
400979	24 G x 3 ½ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.55 x 88.9	Thin	San Agustín	400965	405104
400978	22 G x 3 ½ in. BD [®] Whitacre Spinal NRFit [™] Needle	0.7 x 88.9	Regular	San Agustín	400964	401995

SSU = Sterile, single-use

BNS = Bulk non-sterile



BD® Quincke spinal NRFit™ needles

Code NRFit SSU	Description	Metric dimensions (mm)	Wall	BD manufacturing site	Analogous BNS code	Analogous Luer code
BD® Quincke Spinal NRFit™ Needles						
400932	27 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	0.4 x 88.9	Regular	San Agustín	400988	405259
400931	26 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	0.45 x 88.9	Regular	San Agustín	400091	405258
400063	25 G x 3 in. BD® Quincke Spinal NRFit™ Needle	0.5 x 76.2	Regular	San Agustín	-	-
400089	25 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	0.5 x 88.9	Regular	San Agustín	400080	405257
400219	25 G x 4 1/16 in. BD® Quincke Spinal NRFit™ Needle	0.5 x 119.1	Regular	Juncos	-	-
400930	23 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	0.6 x 88.9	Thin	San Agustín	-	405240
400072	22 G x 1 ½ in. BD® Quincke Spinal NRFit™ Needle	0.7 x 38.1	Regular	San Agustín	400071	405254
400073	22 G x 2 ½ in. BD® Quincke Spinal NRFit™ Needle	0.7 x 63.5	Regular	San Agustín	400064	405244
400077	22 G x 3 in. BD® Quincke Spinal NRFit™ Needle	0.7 x 76.2	Regular	San Agustín	-	405255
400079	22 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	0.7 x 88.9	Regular	San Agustín	400074	405256
400217	22 G x 5 in. BD® Quincke Spinal NRFit™ Needle	0.7 x 127	Regular	Juncos	400075	400366
400218	22 G x 7 in. BD® Quincke Spinal NRFit™ Needle	0.7 x 177.8	Regular	Juncos	400068	400367
400927	20 G x 1 ½ in. BD® Quincke Spinal NRFit™ Needle	0.9 x 38.1	Regular	San Agustín	-	-
400928	20 G x 3 in. BD® Quincke Spinal NRFit™ Needle	0.9 x 76.2	Regular	San Agustín	-	405252
400929	20 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	0.9 x 88.9	Regular	San Agustín	400084	405253
400216	20 G x 6 in. BD® Quincke Spinal NRFit™ Needle	0.9 x 152.4	Regular	Juncos	-	400368
400925	19 G x 3 in. BD® Quincke Spinal NRFit™ Needle	1.1 x 76.2	Regular	San Agustín	-	-
400926	19 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	1.1 x 88.9	Regular	San Agustín	-	405250
400065	18 G x 3 in. BD® Quincke Spinal NRFit™ Needle	1.2 x 76.2	Regular	San Agustín	-	-
400070	18 G x 3 ½ in. BD® Quincke Spinal NRFit™ Needle	1.2 x 88.9	Regular	San Agustín	400069	405248
400215	18 G x 6 in. BD® Quincke Spinal NRFit™ Needle	1.2 x 152.4	Regular	Juncos	-	400370

SSU = Sterile, single-use
BNS = Bulk non-sterile



BD[®] spinal NRFit[™] sets

BD[®] Quincke and Whitacre spinal NRFit[™] needle set

Code NRFit SSU	Description	Metric dimensions (mm)	Wall	BD manufacturing site	Analogous Luer code
BD [®] Spinal NRFit [™] Needle Sets					
400227	27 G x 3 ½ in. BD [®] Whitacre Spinal NRFit [™] Needle Set with 22 G BD [®] Spinal Introducer NRFit [™] Needle	0.4 x 88.9 and 0.7 x 31.8	Regular	San Agustín	405075
400230	27 G x 4 ¼ in. BD [®] Whitacre Spinal NRFit [™] Needle Set with 22 G BD [®] Spinal Introducer NRFit [™] Needle	0.4 x 103.2 and 0.7 x 31.8	Regular	San Agustín	405113
400228	25 G x 3 ½ in. BD [®] Whitacre Spinal NRFit [™] Needle Set with 20 G BD [®] Spinal Introducer NRFit [™] Needle	0.5 x 88.9 and 0.9 x 31.8	Thin	San Agustín	405076
400229	25 G x 4 ¼ in. BD [®] Whitacre Spinal NRFit [™] Needle Set with 20 G BD [®] Spinal Introducer NRFit [™] Needle	0.5 x 103.2 and 0.9 x 31.8	Thin	San Agustín	405112
400226	27 G x 3 ½ in. BD [®] Quincke Spinal NRFit [™] Needle Set with 22 G BD [®] Spinal Introducer NRFit [™] Needle	0.4 x 88.9 and 0.7 x 31.8	Regular	San Agustín	405069
400225	26 G x 3 ½ in. BD [®] Quincke Spinal NRFit [™] Needle Set with 20 G BD [®] Spinal Introducer NRFit [™] Needle	0.45 x 88.9 and 0.9 x 31.8	Regular	San Agustín	405065
400224	25 G x 3 ½ in. BD [®] Quincke Spinal NRFit [™] Needle Set with 20 G BD [®] Spinal Introducer NRFit [™] Needle	0.5 x 88.9 and 0.9 x 31.8	Regular	San Agustín	405084

BD[®] syringes NRFit[™] Lok and BD[®] syringes NRFit[™] Slip

The BD[®] Syringe NRFit[™] Lok and BD[®] Syringe NRFit[™] Slip are intended for neuraxial use by healthcare professionals for aspiration/injection of fluids.

BD[®] syringes NRFit[™] Lok and Slip

Code NRFit SSU	Description	Scale markings		Dead space (maximum, without needle)	Analogous BNS code	Analogous Luer code
		Major graduations	Minor graduations			
BD [®] Syringes NRFit [™]						
400050	3 mL BD [®] Syringe NRFit [™] Lok	0.5 mL	0.1 mL	0.07 mL	400048	309658
400172	3 mL BD [®] Syringe NRFit [™] Slip	0.5 mL	0.1 mL	0.07 mL	400166	-
400173	5 mL BD [®] Syringe NRFit [™] Lok	1 mL	0.2 mL	0.075 mL	400167	309649
400051	5 mL BD [®] Syringe NRFit [™] Slip	1 mL	0.2 mL	0.075 mL	400049	-
400174	10 mL BD [®] Syringe NRFit [™] Lok	1 mL	0.2 mL	0.075 mL	400168	305959
400175	10 mL BD [®] Syringe NRFit [™] Slip	1 mL	0.2 mL	0.075 mL	400169	302146
400182	20 mL BD [®] Syringe NRFit [™] Lok	5 mL	1 mL	0.15 mL	400178	300629
400183	50 mL BD [®] Syringe NRFit [™] Lok	10 mL	1 mL	0.20 mL	400179	300865

BD[®] Blunt needles NRFit[™]

The BD[®] Blunt Fill NRFit[™] and the BD[®] Blunt Filter NRFit[™] Needles are used for aspiration of fluids from vials and ampoules.

BD[®] Blunt needles NRFit[™] Fill and Filter

Code NRFit SSU	Description	Needle dimensions	Analogous Luer code
BD [®] Blunt needles NRFit [™] Fill and Filter			
400066	BD [®] Blunt Filter NRFit [™] Needle	1.2 mm x 40 mm (18 G x 1 ½ in.)	305211
400067	BD [®] Blunt Fill NRFit [™] Needle	1.2 mm x 40 mm (18 G x 1 ½ in.)	303129

SSU = Sterile, single-use

BNS = Bulk non-sterile

BD Switzerland Sàrl, Terre Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland

bd.com

BD, the BD logo and Luer-Lok are trademarks of Becton, Dickinson and Company or its affiliates.
NRFit is a trademark of GEDSA, used with permission. © 2023 BD. All rights reserved. BD-82946 (02/23)

