



Animal Health Catalogue

IV Lines Flushing



BD PosiFlush™ Prefilled Flush Syringe [NaCl 0.9 %]

- Ready-to-use, prefilled flushing device containing sterile sodium chloride (isotonic NaCl 0.9 %) specially designed for flushing in-situ vascular access devices
- Special design of syringe barrel. The 3ml, 5ml and 10ml syringes all have a diameter of a 10ml syringe to ensure low pressure for reduced risk of catheter or vein damage during the flushing process
- BD Luer-Lok™ tip for secure connection to the catheter
- Designed to reduce syringe-induced blood reflux, to help maintain catheter patency
- Size and shape of the tip cap allow easy handling, and help avoid touch contamination of the syringe luer connection
- Eliminates the use of needles



BD PosiFlush™ XS (Externally Sterile) - for use in critical care areas for sterile procedures and within a sterile field

Description	Shelf Box / Shipping Case	Reference
Posiflush XS 3ml [NaCl 0.9%]	30 / 240	306570
Posiflush XS 5ml [NaCl 0.9%]	30 / 240	306571
Posiflush XS 10ml [NaCl 0.9%]	30 / 240	306572

BD PosiFlush™ SP - Sterile Fluid Path, for general IV flush procedures

Description	Shelf Box / Shipping Case	Reference
Posiflush SP 3ml [NaCl 0.9%]	30 / 480	306573
Posiflush SP 5ml [NaCl 0.9%]	30 / 480	306574
Posiflush SP 10ml [NaCl 0.9%]	30 / 480	306575

BD Switzerland Sàrl, Terre Bonne Park – A4 Route de Crassier 17, 1262 Eysins, Switzerland
Tel: +41 21 556 30 00

[bd.com](https://www.bd.com)

