

BD BACTEC™ Synapsys™ Adaptor

September 2020

Critical and/or Security Patches June 2020.

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related through June 2020.

The patches include updates that may adversely affect BD BACTEC Synapsys Adaptor. These patches will not be and should not be installed on BD BACTEC Synapsys Adaptor at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Patch Description	Patch ID	Month
------------	-------------------	----------	-------

Notes



BD BACTEC™ Synapsys™ Adaptor

September 2020

Critical and/or Security Patches for June 2020

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related through June 2020. These patches were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Service Stack Update for Windows 10 IoT for x64-based Systems	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates.	4562561	June 2020
Security Update for Adobe Flash Player	This security update resolves vulnerabilities in Adobe Flash Player that is installed on any of the Windows operating systems	4561600	June 2020
Cumulative Security Update for Windows 10 IoT for x64-based Systems	<p>This security update includes quality improvements. Key changes include:</p> <ul style="list-style-type: none">• Addresses an issue that prevents users from updating .msi files from a network folder.• Addresses an issue that prevents you from installing some .msi apps. This occurs when a device is managed by a Group Policy that redirects the AppData folder to a network folder.• Security updates to the Microsoft Scripting Engine, Internet Explorer, Windows App Platform and Frameworks, Windows Media, Windows Kernel, Microsoft Graphics Component, Windows Input and Composition, Windows Shell, the Microsoft Store, Windows Cloud Infrastructure, Windows Fundamentals, Windows Authentication, Windows Cryptography, Windows Virtualization, Windows Storage and Filesystems, Windows File Server and Clustering, Windows Hybrid Storage Services, the	4561616	June 2020



Advancing the world of health

	Windows Update Stack, and the Microsoft JET Database Engine.		
--	--	--	--

Notes

1. Service Stack update 4562561 must be installed prior to other updates listed.