

BD Assurity Linc

Date of Critical or Security Patches: May 2019

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for May 2019. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Notes
Intel microcode updates	This security update resolves various vulnerabilities in Windows Server 2016, Windows 10, version 1607.	<a href="#">KB4346087</a>	1
Servicing stack update for Windows Server 2016 and Windows 10, version 1607: February 12, 2019	This security This update makes quality improvements to the servicing stack component that installs Windows updates.	<a href="#">KB4485447</a>	1

### Notes

1. Tested against BD Assurity Linc configured with Windows 10 IoT Professional 64-bit build 1607



Advancing the world of health