

BD Assurity Linc™

December 2020

Critical and/or Security Patches December 2020.

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related through September 2020.

The patches include updates that may adversely affect BD Assurity Linc. These patches will not be and should not be installed on BD Assurity Linc at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Patch Description	Patch ID	Month
------------	-------------------	----------	-------

Notes



BD Assurity Linc™

December 2020

Critical and/or Security Patches December 2020

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related through September 2020. These patches were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Service Stack Update for Windows 10 IoT for x64-based Systems	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates.	4576750	September 2020
Cumulative Security Update for Windows 10 IoT for x64-based Systems	This security update includes quality improvements. Key changes include: <ul style="list-style-type: none">• Provides the ability to set a Group Policy that displays only the domain and username when you sign in.• Updates time zone information for the Yukon, Canada.• Addresses an issue that prevents the Event Viewer from properly saving filtered events.• Addresses an issue with object performance counters.• Addresses an issue that sometimes prevents AppLocker from running an application whose publisher rule allows it to run.• Addresses an issue that causes delays during shutdown when running the Microsoft Keyboard Filter Service.• Addresses a security vulnerability issue with user proxies and HTTP-based intranet servers. After installing this update, HTTP-based intranet servers cannot leverage a user proxy by default to detect updates. Scans using these	4577015	September 2020



	<p>servers will fail if the clients do not have a configured system proxy. If you must leverage a user proxy, you must configure the behavior using the Windows Update policy "Allow user proxy to be used as a fallback if detection using system proxy fails." This change does not affect customers who secure their Windows Server Update Services (WSUS) servers with the Transport Layer Security (TLS) or Secure Sockets Layer (SSL) protocols. For more information, see Improving security for devices receiving updates via WSUS.</p> <ul style="list-style-type: none"> • Addresses an issue with evaluating the compatibility status of the Windows ecosystem to help ensure application and device compatibility for all updates to Windows. • Addresses an issue that might prevent a client from reconnecting to a previous RemoteApp session unless an administrator closes the session on the server. • Security updates to Windows App Platform and Frameworks, Microsoft Graphics Component, Windows Input and Composition, Windows Media, Windows Shell, the Microsoft Store, Windows Cloud Infrastructure, Windows Fundamentals, Windows Kernel, Windows Virtualization, Windows Storage and Filesystems, the Microsoft Scripting Engine, and the Microsoft JET Database Engine. 		
--	--	--	--

Notes

1. Service Stack update 4576750 must be installed prior to other updates listed.

