

# Security Patches: BD Assurity Linc™

September 2023

BD has identified patches from Microsoft® that have been identified as critical or security related through July 2023 for the BD Assurity Linc™ on Windows LTSC 1809. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

## Microsoft® patches

Patch name	Description	Patch ID	Notes
Security update for Windows 10 1809.	This update addresses security issues in the Windows operating system.	<a href="#">KB5028168</a>	
Cumulative security update for .NET framework 3.5 and 4.7.2 for Win 10 1809.	This security update addresses multiple remote code execution, elevation of privilege and denial of service vulnerabilities.	<a href="#">KB5027131</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

