

BD Assurity Linc™

Date of Critical or Security Patches: September 2021

Critical Security Patches: April 2021-July 2021

Microsoft® & Third-Party Patches for Windows 10 1607

BD has identified patches from Microsoft® that have been identified as critical or important security related updates thru July 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Notes
Cumulative Security Update	<ul style="list-style-type: none">▪ Updates to remove Adobe Flash from your device.▪ Updates for verifying usernames and passwords.▪ Updates to improve security when Windows performs basic operations.	KB5004238	
Servicing Stack Update	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates.	KB5001402	1

Notes:

Servicing Stack update KB5001402 must be installed prior to other updates.



Critical Security Patches: February 2021-July 2021

Microsoft® & Third-Party Patches for Windows 10 1809

BD has identified patches from Microsoft® that have been identified as critical or important security related updates thru July 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Notes
Cumulative Security Update for .NET Framework	This security update addresses a denial of service vulnerability in .NET Framework	KB4601060	
Cumulative Security Update	<ul style="list-style-type: none">Updates to improve security when Windows performs basic operations.Updates for verifying usernames and passwords.	KB5004244	
Servicing Stack Update	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates.	KB5003711	1

Notes:

Servicing Stack update KB5003711 must be installed prior to other updates.