

Security Patches: BD Assurity Linc™

July 2025

BD has identified patches from Microsoft® that have been identified as critical or security related through May 2025 for BD Assurity Linc™ on Windows LTSC 1809. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Microsoft® patches

Patch name	Description	Patch ID	Notes
Cumulative Update for .NET Framework 3.5 and 4.7.2 for Windows 10, version 1809	This update is cumulative and contains all previously released security improvements.	KB5054695	
Cumulative security Update for Windows 10 Version 1809	This update addresses security issues in the Windows operating system	KB5058392	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

