

Security Patches: BD BACTEC™ FX

January 2025

BD has identified patches from Microsoft® that have been identified as critical or security related through October 2024. An update to the anti-malware software provisioned in BD BACTEC FX was also validated recently and can be applied by BD trained service engineers.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

Microsoft® patches

Patch name	Description	Patch ID	Notes
Servicing Stack Update	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5043124	Must be installed before other updates.
Windows 10 1607 Cumulative Security update.	Addresses several security issues for the Windows operating system.	KB5044293	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

