

Security Patches: BD BACTEC™ FX40

September 2025

BD has identified patches from Microsoft® that have been identified as critical or security related through July 2025. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

Microsoft® patches

| Patch name | Description | Patch ID | Notes |
|---|---|---------------------------|---|
| Servicing Stack Update | This update makes quality improvement to the servicing stack, which is the component that installs Windows updates. | KB5062799 | Must be installed before other updates. |
| Windows 10 1607 Cumulative Security update. | Addresses several security issues for the Windows operating system. | KB5063871 | |

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

