

BD Product Name: **BD Care Coordination Engine**

Date of Critical or Security Patches: December 2020

Abstract: Critical or Security Patches – December 2020

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for December 2020.

The patches include updates that may adversely affect BD Care Coordination Engine. These patches will not be and should not be installed on BD Care Coordination Engine at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
Security Update for Microsoft Silverlight	This security update to Silverlight includes fixes outlined in KB 4013867. This update is backward compatible with web applications built using previous versions of Silverlight.	KB4013867	None
2020-02 Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article.	KB4537810	None
2020-02 Security Only Quality Update for Windows Server 2008 for x64-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can	KB4537822	None



	help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update		
2020-05 Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 for x64	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update	KB4556402	None
2020-10 Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 for x64	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update	KB4566469	None