

BD Product Name: **Care Coordination Engine**

Date of Critical or Security Patches: January 2021

Abstract: Critical or Security Patches – January 2021

### **Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for January 2021.

The patches include updates that may adversely affect BD CCE. These patches should be installed CCE Windows Server 2008 standard.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

<b>Patch Name</b>	<b>Description</b>	<b>Patch ID</b>	<b>Notes</b>
This security update to Silverlight includes fixes outlined in Security update for Server 2008 Standard	This security update to Silverlight includes fixes outlined in Security update KB4013867. Non-security related changes are also included in this release and have been documented in the Silverlight Release History. This update is backward compatible with web applications built using previous versions of Silverlight.	<a href="#">KB4013867</a>	1
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems for 2008 Server	This security update includes improvements and fixes that were a part of update KB4534303 (released January 14, 2020) and addresses the following issues: Security updates to Microsoft Graphics Component, Windows Input and Composition, Windows Shell, Windows Fundamentals, Windows Cryptography, Windows Hyper-	<a href="#">KB4537810</a>	1



	V, Windows Core Networking, Windows Peripherals, Windows Network Security and Containers, Windows Storage and Filesystems, the Microsoft Scripting Engine, and Windows Server.		
This security update includes quality improvements for 2008 Server.	Security updates to Microsoft Graphics Component, Windows Input and Composition, Windows Shell, Windows Fundamentals, Windows Cryptography, Windows Hyper-V, Windows Core Networking, Windows Peripherals, Windows Network Security and Containers, Windows Storage and Filesystems, and Windows Server.	<a href="#">KB44537822</a>	1
Microsoft .NET Framework 2.0 Microsoft .NET Framework 3.0 Microsoft .NET Framework 4.5.2 Microsoft .NET Framework 4.6 for 2008 Server.	Security updates to Microsoft Graphics Component, Windows Input and Composition, Windows Shell, Windows Fundamentals, Windows Cryptography, Windows Hyper-V, Windows Core Networking, Windows Peripherals, Windows Network Security and Containers, Windows Storage and Filesystems, and Windows Server.	<a href="#">KB4566469</a>	1
Microsoft .NET Framework 2.0 Microsoft .NET Framework 3.0 Microsoft .NET Framework 4.5.2 Microsoft .NET Framework 4.6 for 2008 Server.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB4556406</a>	1

Microsoft .NET Framework 2.0 Microsoft .NET Framework 3.0 Microsoft .NET Framework 4.5.2 Microsoft .NET Framework 4.6 for 2008 Server.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB4575904</a>	1
Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB4570503</a>	1
Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB4579980</a>	1

<p>Security Only Update for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008</p>	<p>A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.</p>	<p><a href="#">KB4580470</a></p>	<p>1</p>
---	---	----------------------------------	----------

1. Note 1: These available updates are not able to install due to an expired license issue on windows 2008 Server 64-bit operating system.