

Security Patches: Care Coordination Engine

June 2022

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2022 will not install on CCE.

The patch includes updates that may adversely affect CCE. This patch will not be and should not be installed on CCE at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

Microsoft® patches

Patch name	Description	Patch ID	Notes
2022-06 Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5014752	1
2022-06 Security Monthly Quality Rollup for Windows Server 2008 R2 for x64-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5014748	1
2022-06 Security Only Quality Update for Windows Server 2008 for x64-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5014743	1

2022-06

Security Only
Quality Update
for Windows
Server 2008 R2
for x64-based
Systems

A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.

[KB5014742](#)

1

Note:

1. These available updates are not able to install due to an expired license issue on DoD 1.4.1, CCER 1.3.1, CCE 1.4.1 Windows 2008 Server 64-bit operating system. Hence, we were unable to test these.