

Security Patches: Care Coordination Engine

February 2024

BD has identified patches from Microsoft® that have been identified as critical or security related for February 2024. The patches include updates that may adversely affect BD CCE. These patches should not be installed on CCE at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Microsoft® patches

Patch name	Description	Patch ID	Notes
Microsoft Edge-Extended Stable Channel Version 122 Update for x64 based Editions (Build 122.0.2365.66)	This update provides the latest feature and quality updates to Microsoft Edge. These updates are designed to be deployed with Windows Server Update Services and Configuration Manager. See the "More Information" link for a list of supported operating systems.	122.0.2365.66	1
2024-02 Security and Quality Rollup for .NET Framework 3.5, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8 for Windows Server 2012 R2 for x64	There are no new security improvements in this release. This update is cumulative and contains all previously released security improvements.	KB5034689	1
2024-02 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems	This cumulative security update includes improvements that are part of update KB5034171.	KB5034819	1

2024-02 Servicing Stack
Update for Windows Server
2012 R2 for x64-based
Systems

This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. Servicing stack updates (SSU) makes sure that you have a robust and reliable servicing stack so that your devices can receive and install Microsoft updates.

[KB5034866](#)

1

Note: -

1. Patches (2012R2) could not be validated on DOD server due to environment issues.

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

[bd.com](#)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

