

Security Patches: Care Coordination Engine

June 2023

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2023. The patches include updates that may adversely affect BD CCE. These patches should not be installed on CCE at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Microsoft® patches			
Patch name	Description	Patch ID	Notes
Microsoft Edge-Stable Channel Version 114 Update for x64 based Editions (Build 114.0.1823.43)	This update provides the latest feature and quality updates to Microsoft Edge. These updates are designed to be deployed with Windows Server Update Services and Configuration Manager. See the "More Information" link for a list of supported operating systems.	Microsoft Edge-Stable Channel Version 114 Update for x64 based Editions (Build 114.0.1823.43)	1

Note: -

1. Patch could not be verified due to system environment issues.

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

