

Security Patches: Care Coordination Engine

October 2023

BD has identified patches from Microsoft® that have been identified as critical or security related for October 2023. The patches include updates that may adversely affect BD CCE. These patches should not be installed on CCE at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Microsoft® patches			
Patch name	Description	Patch ID	Notes
2023-10 Cumulative Update for Windows Server 2016 for x64-based Systems (KB5031362)	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5031362	1

Note: -

1. Patch(es) could not be verified due to system application issues.

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

[bd.com](#)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

