

BD Product Name: **Care Coordination Engine**

Date of Critical or Security Patches: Nov 2021

Abstract: Critical or Security Patches – Nov 2021

### **Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for Nov 2021.

The patches include updates that may adversely affect BD CCE. These patches should not be installed CCE Windows Server 2008 standard.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

<b>Patch Name</b>	<b>Description</b>	<b>Patch ID</b>	<b>Notes</b>
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems.	The Microsoft Security Response Center (MSRC) investigates all reports of security vulnerabilities affecting Microsoft products and services, and provides the information here as part of the ongoing effort to help you manage security risks and help keep your systems protected.	<a href="#">KB5007263</a>	1
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems.	Addresses a known issue that might prevent devices from downloading and installing printer drivers when the devices attempt to connect to a network printer for the first time. We have observed this issue on devices that access printers using a print server that uses HTTP connections.	<a href="#">KB5007246</a>	1
Security and Quality Rollup for .NET Framework 4.5.2 for Windows Server 2008 SP2 for x64.	Install this update to resolve issues in Windows. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article for more information. After you install this item, you may have to restart your computer..	<a href="#">KB5007167</a>	1



1. Note 1: This available update is not able to install due to an expired license issue on DoD 1.4.1, CCER 1.3.1, CCE 1.4.1 Windows 2008 Server 64-bit operating system.

