

BD Product Name: **Care Coordination Engine**

Date of Critical or Security Patches: Sept 2021

Abstract: Critical or Security Patches – Sept 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for Sept 2021.

The patches include updates that may adversely affect BD CCE. These patches should not be installed CCE Windows Server 2008 standard.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Description	Patch ID	Notes
Windows Malicious Software Removal Tool x64.	After the download, this tool runs one time to check your computer for infection by specific, prevalent malicious software (including Blaster, Sasser, and Mydoom) and helps remove any infection that is found. If an infection is found, the tool will display a status report the next time that you start your computer. A new version of the tool will be offered every month. If you want to manually run the tool on your computer, you can download a copy from the Microsoft Download Center, or you can run an online version from microsoft.com. This tool is not a replacement for an antivirus product. To help protect your computer, you should use an antivirus product.	KB890830	1
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can	KB5004955	1



	help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.		
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system..	KB4537810	1
Cumulative Security Update for Internet Explorer 9 for Windows Server 2008 for x64-based systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5005563	1
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you	KB5005606	1

	may have to restart your system.		
Security Only Quality Update for Windows Server 2008 for x64-based Systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5005618	1

1. Note 1: This available update is not able to install due to an expired license issue on DoD Windows 2008 Server 64-bit operating system.